


Compact spectrometers with
built-in Hamamatsu image
sensor, optical element, etc.

Mini- spectrometers



Related product



FTIR engines (FT-NIR spectrometer)

Portable NIR
spectroscopic modules



MEMS-FPI spectroscopic modules

Ultra-compact near infrared
spectrum sensor with
MEMS-FPI tunable filter

Mini-spectrometers

[Home](#)[What are mini-spectrometers?](#)[Mini-spectrometers lineup](#)[Spectrometer heads](#)[Technologies](#)[Application examples](#)[Technical note](#)[Accessories](#)

We have more than 20 different mini-spectrometers for the ultraviolet to near infrared regions.

What are mini-spectrometers?


Mini-spectrometers are small spectrometers (polychromators) with an integrated optical system, image sensor, and driver circuit. They are portable devices that make them possible to do real-time measurement on-site.



Applications

- Color measurement
- Sugar content measurement
- Film thickness measurement
- Plastic screening
- Fluorescence measurement
- Environmental analysis
- Mobile measuring devices

● Example of mini-spectrometer optical system



KACCC0256EA




Ultra-small spectrometer heads (without a driver circuit) are also available.

Series	Products	Spectral response range (nm)															
		UV		Visible				Near infrared									
		200	400	600	800	1000	1200	1400	1600	1800	2000	2200	2400	2600			
For ultraviolet range	High sensitivity C10082CA High resolution C10082CAH		200 to 800														
	Wide dynamic range C10082MD		200 to 800														
	High resolution C16449MA-02		200 to 600														
	High sensitivity C9404CA High resolution C9404CAH		200 to 400														
For ultraviolet to near infrared range	High sensitivity C16449MA-01		190 to 1100														
For visible range	High sensitivity C10083CA High resolution C10083CAH		320 to 1000														
	Wide dynamic range C10083MD High sensitivity C11697MB		320 to 1000														
	High sensitivity C13555MA		340 to 830														
For visible to near infrared range	High near IR sensitivity C9405CC High sensitivity C13053MA		500 to 1100														
For near infrared range	Non-cooled type C11482GA Cooled type C9913GC		900 to 1700														
	Cooled type C9914GB		1100 to 2200														
	Cooled type C11118GA		900 to 2550														
	Compact type C14486GA		950 to 1700														
For Raman spectroscopy	High resolution C14214MA		790 to 1050														

Note: See [P.13](#) for details on spectrometer heads.


Spectral response



* A/D count when constant light level enters optical fiber
(Fiber core diameter: 600 μm , assuming no attenuation in optical fiber)

KACCB0137EQ







Spectral resolution vs. wavelength (typical example)



KACCB0139EQ


For ultraviolet range

These are products with sensitivity in the ultraviolet range.

Type no.	Type	Spectral response range (nm)							Spectral resolution typ. (nm)	S/N max.	External power supply	Built-in image sensor	Size (mm)	Photo
		UV		Visible			Near infrared							
		200	400	600	800									
C10082CA	High sensitivity				200 to 800				4	446 : 1	+5 V	Back-thinned CCD S10420-1106-01	95 × 92 × 76	
C10082CAH	High resolution				200 to 800				1	446 : 1	+5 V	Back-thinned CCD S10420-1106-01	95 × 92 × 76	
C10082MD	Wide dynamic range				200 to 800				4	4390 : 1	Not required (USB bus power only)	CMOS linear image sensor S8378-1024Q	94 × 90 × 55	
C9404CA	High sensitivity		200 to 400						2	446 : 1	+5 V	Back-thinned CCD S10420-1006-01	125.7 × 115.7 × 75	
C9404CAH	High resolution		200 to 400						1	446 : 1	+5 V	Back-thinned CCD S10420-1006-01	125.7 × 115.7 × 75	
C16449MA-02	High sensitivity		200 to 600						0.45	316 : 1	Not required (USB bus power only)	High-sensitivity CMOS image sensor	80 × 75 × 25	






For ultraviolet to near infrared range

It is a wide spectral response type with sensitivity extending from the UV region.

Type no.	Type	Spectral response range (nm)										Spectral resolution typ. (nm)	S/N max.	External power supply	Built-in image sensor	Size (mm)	Photo
		UV		Visible				Near infrared									
		200	400	600	800	1000	1200										
C16449MA-01	High resolution											1	316 : 1	Not required (USB bus power only)	High-sensitivity CMOS image sensor	80 × 75 × 25	



For visible range

These are products with sensitivity in the visible range.

Type no.	Type	Spectral response range (nm)								Spectral resolution typ. (nm)	S/N max.	External power supply	Built-in image sensor	Size (mm)	Photo
		UV		Visible			Near infrared								
		200	400	600	800	1000									
C10083CA	High sensitivity			320 to 1000						5	446 : 1	+5 V	Back-thinned CCD S10420-1106-01	95 × 92 × 76	
C10083CAH	High resolution			320 to 1000						1	446 : 1	+5 V	Back-thinned CCD S10420-1106-01	95 × 92 × 76	
C10083MD	Wide dynamic range			320 to 1000						5	4390 : 1	Not required (USB bus power only)	CMOS linear image sensor S8378-1024Q	94 × 90 × 55	
C11697MB	High sensitivity			320 to 1000						5	260 : 1	Not required (USB bus power only)	High-sensitivity CMOS linear image sensor S11639	94 × 90 × 55	
C13555MA	High sensitivity			340 to 830						2.3	230 : 1	Not required (USB bus power only)	High-sensitivity CMOS linear image sensor	80 × 60 × 12	






For visible to near infrared range

These are products with sensitivity in the visible to near infrared range.

Type no.	Type	Spectral response range (nm)											Spectral resolution typ. (nm)	S/N max.	External power supply	Built-in image sensor	Size (mm)	Photo
		UV		Visible			Near infrared											
		200	400	600	800	1000	1200											
C9405CC	High near IR sensitivity												4	446 : 1	+5 V	Back-thinned CCD S16010-1006	125.7 × 115.7 × 75	
C13053MA	High sensitivity												2.5	230 : 1	Not required (USB bus power only)	High-sensitivity CMOS linear image sensor	80 × 60 × 12	


For near infrared range

These are products with sensitivity in the near infrared range.

Type no.	Type	Spectral response range (nm)												Spectral resolution typ. (nm)	S/N max.	External power supply	Built-in image sensor	Size (mm)	Photo
		Near infrared																	
		800	1000	1200	1400	1600	1800	2000	2200	2400	2600								
C11482GA	Non-cooled type													5	7700 : 1	Not required (USB bus power only)	InGaAs linear image sensor G9204-512DA	38.5 × 106 × 86	
C9913GC	Cooled type													5	6100 : 1	+5 V, +12 V	InGaAs linear image sensor G9204-512SA	142 × 218 × 82	
C9914GB	Cooled type													6	6100 : 1	+5 V, +12 V	InGaAs linear image sensor	142 × 218 × 82	
C11118GA	Cooled type													15	7700 : 1	+5 V, +12 V	InGaAs linear image sensor G9208-256WB-02	142 × 218 × 82	
C14486GA	Compact type													5	6900 : 1	Not required (USB bus power only)	InGaAs linear image sensor	80 × 60 × 12	





For Raman spectroscopy

It is a product with sensitivity in the Raman spectroscopy.

Type no.	Type	Spectral response range (nm)								Spectral resolution typ. (nm)	S/N max.	External power supply	Built-in image sensor	Size (mm)	Photo
		UV		Visible			Near infrared								
		200	400	600	800	1000									
C14214MA	High resolution								790 to 1050	0.4	230 : 1	Not required (USB bus power only)	High-sensitivity CMOS linear image sensor	100 × 60 × 12	

Spectrometer heads


The compact spectrometer heads (without a driver circuit) integrate an optical system and an image sensor.

Type no.	Type	Spectral response range (nm)								Spectral resolution typ. (nm)	S/N max.	Built-in image sensor	Size (mm)	Photo
		UV		Visible			Near infrared							
		200	400	600	800	1000								
C16767MA	For ultraviolet range	190 to 440								5.5	293 : 1	High-sensitivity CMOS linear image sensor	20.1 × 12.5 × 10.1	
C12666MA	Wide dynamic range	340 to 780								12	5300 : 1	CMOS linear image sensor	20.1 × 12.5 × 10.1	
C12880MA	High sensitivity	340 to 850								12	291 : 1	High-sensitivity CMOS linear image sensor	20.1 × 12.5 × 10.1	
C11708MA	For near IR	640 to 1050								15	5300 : 1	CMOS linear image sensor	27.6 × 16.8 × 13	

Optical system in the compact spectrometer heads

The miniaturization of the spectrometer head has been achieved by employing a CMOS image sensor with a slit formed by etching and a grating fabricated by nanoimprinting within the optical system.


C12666MA, C12880MA, C16767MA



KACCC1035EB

The metal package provides high humidity resistance (C12666MA, C12880MA). Low cost is achieved because it is a hollow type.

C11708MA



KACCC0922EC


The glass used does not expand easily with rising temperatures, so the temperature dependency of the wavelength is extremely small.

Mini-spectrometers, employ MOEMS (micro-opto-electro-mechanical-systems) technology, combining an image sensor / optical system and MEMS.


MOEMS technologies

Image sensors

- Uses one of Hamamatsu image sensor lineup to support various wavelengths
- Available with custom design



▲ CCD image sensor



▲ High-sensitivity CMOS linear image sensor




▲ TE-cooled InGaAs linear image sensor




▲ IR-enhanced CMOS linear image sensor

Optical system


- Optical design suitable for spectrometers
- Optical simulation



MEMS



▲ Grating that uses nanoimprint




▲ Image sensor with a through-hole slit



Software


Supports various
communication interfaces
(e.g., USB)

Evaluation software
available ►



Circuits

- Unique driver circuits
- Evaluation circuits available
for spectrometer heads




Mini-spectrometers

Application examples

Mini-spectrometers can be incorporated into a variety of devices and are used in a wide range of applications.

Color measurement (e.g., LED light source)



KACCC0796EA

A mini-spectrometer is used to perform spectral measurement and inspect LEDs or the like.


Sugar content measurement



KACCC0797EA

A mini-spectrometer is used in applications such as handy brix meters, which measure sugar content by absorbance.


Plastic screening



KACCC0601EB

Plastic screening is performed by using the fact that when near infrared light is directed at plastic, the wavelengths that are absorbed varies depending on the material.


Display color measurement



KACCC0599EC

The emission spectrum of LCDs is monitored with a mini-spectrometer.


Film thickness measurement



KACCC0600EB

White light interferometry is used to measure the spectrum peak count, film refractive index, and film thickness from the light incident angle.

Environmental analysis



KACCC0798EB

A mini-spectrometer is used in environmental analysis of water, soil, and the like.

Mini-spectrometers

Accessories

Accessories for mini-spectrometers (sold separately) are available.

Input optical fibers A16962 series, A16963 series

UV/visible optical fiber (UV resistant) and visible/NIR optical fiber are available.

Type no.	Product name	Core diameter (μm)	Minimum bend radius (mm)	Specification
A16962-01	Ultraviolet/visible optical fiber (UV resistant)	600	132	NA=0.22 1.5 m in length, with SMA905D connector on each end Operating temperature: 0 °C to +60 °C Storage temperature: -10 °C to +70 °C
A16962-02		800	176	
A16963-01	Visible/near infrared optical fiber	600	132	
A16963-02		800	176	

External trigger coaxial cables A10670, A12763

Cable	Applicable mini-spectrometers	Length (m)
A10670	C9404CA, C9404CAH, C10082CA, C10082CAH, C10082MD, C10083CA, C10083CAH, C10083MD, C11118GA, C11697MB, C11482GA	1.5
A12763	C13555MA, C13053MA, C14486GA, C14214MA, C16449MA-01, C16449MA-02	

2 W xenon flash lamp modules L13651 series



These lamp modules integrate a 2 W xenon flash lamp with a power supply and trigger socket, and are designed to extract maximum performance from the lamp.

Features

- Compact: 42 mm × 42 mm × 37 mm
- Operates on 5 V mobile battery
- Long life: 1×10^9 flash
- Repetition rate: 1250 Hz max.
- Broad spectrum: UV region to middle IR region

[Note: We offer a catalog of xenon flash lamps.](#)

FTIR engines (FT-NIR spectrometer)

Compact FT-NIR spectroscopic modules that can be incorporated into portable measuring instruments and in-line measuring instruments



C1511-01

C16551-01

The Fourier transform infrared spectrometer (FTIR) engines are compact enough to carry in just one hand. A Michelson optical interferometer and a control circuit are built into a palm-sized case. Spectrum and absorbance can be measured by connecting a PC via USB.

Features

- Compact
- High speed 275 frames/s typ. (C16511-01)
- Optical fiber input type
- High S/N
- Suitable for diffuse reflection measurements and absorbance measurements
- Spectral response range: 1100 nm to 2500 nm
- Ethernet compatible (C16511-01)

Applications


- FA, PAT (Process Analytical Technology)
- Bioprocess analysis
- Material analysis
- Farm product and food inspection
- Plastic sorting
- Medicines inspection

Type no.	Spectral response range (nm)													Spectral resolution (nm)
	Near infrared													
	800	1000	1200	1400	1600	1800	2000	2200	2400	2600				
C15511-01 C16511-01														5.7 typ. (λ=1533 nm)

Optical system

The optical interferometer of the FTIR engine consists of a MEMS chip, as well as the light input section, beam splitter, fixed mirror, and photodetector.

Optical system of FTIR engine




MEMS-FPI spectroscopic modules

These compact modules have a built-in MEMS-FPI spectrum sensor and light source.



● Spectral resolution vs. peak transmission wavelength (typical example)







KACCB0764EA

MEMS-FPI spectroscopic modules	Spectral response range (nm)										Spectral resolution (full width at half maximum) max. (nm)
	Near infrared										
	800	1000	1200	1400	1600	1800	2000	2200			
C17552			1350 to 1650								18
C17553				1550 to 1920							21
C17554					1750 to 2150						22

Compact spectrometers for near infrared range

A wide variety of compact spectrometers for the near infrared region are available.

Product name	Type no.	Spectral response range (nm)												Features	Spectroscopic technology	Spectral resolution (nm)	Size (mm)	Photo	
		Near infrared																	
		800	1000	1200	1400	1600	1800	2000	2200	2400	2600								
FTIR engine (FT-NIR spectrometer)	C15511-01														High precision measurement, high wavelength accuracy	MEMS-FTIR	5.7 typ. (λ=1533 nm)	49 × 57 × 76	
FTIR engine (FT-NIR spectrometer)	C16511-01														High precision measurement, high wavelength accuracy, ethernet compatible	MEMS-FTIR	5.7 typ. (λ=1533 nm)	68 × 57 × 76	
MEMS-FPI spectroscopic modules	C17552 C17553 C17554														Compact, suitable for portable devices, excellent high-volume producibility, built-in light source	Fabry-Perot	22 max. (C17554, λ=2150 nm)	74 × 32 × 16	
Mini-spectrometer	C14486GA														High-speed measurement, high sensitivity	Grating	5.0 typ.	80 × 60 × 12	

- [Disclaimer](#)
- [Precautions / Mini- spectrometers](#)

www.hamamatsu.com

- Information described in this material is current as of September 2025.
- Product specifications are subject to change without prior notice due to improvements or other reasons. Before using these products, always contact us for the delivery specification sheet to check the latest specifications.

HAMAMATSU PHOTONICS K.K.

KACC0002E31 Sep. 2025 DN

Optical Semiconductor Sales, HAMAMATSU PHOTONICS K.K.

1126-1 Ichino-cho, Chuo-ku, Hamamatsu City, Shizuoka Pref., 435-8558 Japan, Telephone: (81)53-434-3311, Fax: (81)53-434-5184

U.S.A.: HAMAMATSU CORPORATION: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218

Germany: HAMAMATSU PHOTONICS DEUTSCHLAND GMBH: Arzbergerstr. 10, 82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: HAMAMATSU PHOTONICS FRANCE S.A.R.L.: 19 Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10 E-mail: infos@hamamatsu.fr

United Kingdom: HAMAMATSU PHOTONICS UK LIMITED: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, UK, Telephone: (44)1707-294888, Fax: (44)1707-325777 E-mail: info@hamamatsu.co.uk

North Europe: HAMAMATSU PHOTONICS NORDEN AB: Torshamnsgatan 35, 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se

Italy: HAMAMATSU PHOTONICS ITALIA S.R.L.: Strada della Moia, 1 int. 6 20044 Arese (Milano), Italy, Telephone: (39)02-93 58 17 33, Fax: (39)02-93 58 17 41 E-mail: info@hamamatsu.it

China: HAMAMATSU PHOTONICS (CHINA) CO., LTD.: 1201, Tower B, Jiaming Center, 27 Dongsanhuan Beilu, Chaoyang District, 100020 Beijing, P.R. China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2866 E-mail: hpc@hamamatsu.com.cn

Taiwan: HAMAMATSU PHOTONICS TAIWAN CO., LTD.: 13F-1, No.101, Section 2, Gongdao 5th Road, East Dist., Hsinchu City, 300046, Taiwan(R.O.C) Telephone: (886)3-659-0080, Fax: (886)3-659-0081 E-mail: info@hamamatsu.com.tw

Korea: HAMAMATSU PHOTONICS KOREA CO., LTD.: A-912, 167, Songpa-daero, Seoul, 05855, Korea, Telephone: (82)2-2054-8202, Fax: (82)2-2054-8207 E-mail: sales@hpkcr.co.kr