

HAIGH

Hygienic Solutions



Sluicemaster | **SOLO**[®]
Single occupancy bedpan disposal



www.haigh.co.uk



Hygienic Solutions

Sluicemaster SOLO[®]

Single occupancy bedpan disposal



Eliminating infection transfer across wards

The SOLO[®], installed in the bathroom of each individual hospital isolation room eliminates potential wide area contamination, as only new bedpans enter the room. Current practice with single room bedpan disposal is to use the communal disposer in another multi occupancy ward or treat as clinical waste and double bag. Disposal is then elsewhere which entails moving a used bedpan out of isolation into the general hospital environment.

Besides isolation wards the SOLO[®] can also be used in clinics and other low usage areas*

How it works...

Disposing of just one bedpan at a time, the SOLO[®] is designed for installation in single bed wards. When a used bedpan is ready for disposal the operator walks up to the machine, pushes the lid handle allowing the lid to rise automatically. The bedpan is placed inside, the lid closed and the start button pressed. In under a minute the machine has broken the bedpan up and flushed it away completely.

1. Open lid
2. Insert bedpan
3. Close lid
4. Press button
5. Walk away after washing hands



* SOLO[®] comes with a twelve month guarantee



Sluicemaster SOLO® - Installation & Technical Data

Installation planning

Maximum run to drain 2 metres with a minimum fall of 1:25.

Minimum waste pipe size 38mm (1½").

Use swept bends not 90° degree bends.

Cold water supply

Mains pressure. Ensure supply line is at no point less than 15mm.

Electrical supply

SOLO® is supplied as a single phase machine which must be connected to a 13 Amp fused spur switch (customer supply).

The electrical supply must be to current I.E.C. regulations using the supplied two metre flying lead. The appliance must be connected to a protective earth connection.

Drainage


The machine is fitted with a 38mm (1½") trap inside the cabinet. Ensure the route to the soil stack is shortest with minimum number of bends.

Please refer to the manual supplied with the machine for latest information and installation details as we continuously improve our products over time.

Electrical Supply

Volts	Power	Isolator Fuse	Water Consumption	Cycle Time	Max Loading	Weight (Dry)
240 / 1 / 50	485W	13A	11 ltr / cycle	53 seconds	1 Pulp Product	50Kg

Full dimensions



Notes:

1. Drain outlet terminates at a fitting suitable for a 38mm (1½") dia pipe Stub. Bath trap supplied and fitted internally.
2. Water supply in via 700mm long 10mm bore flexible water inlet pipe with 3/4" bsp female connector supplied.
3. Electrical supply in via 2m long sheathed electrical cable supplied.

All dimensions shown in mm



Hygienic Solutions

www.haigh.co.uk

It is the policy of our company continually to improve our products and accordingly we reserve the right to alter specifications and appearance without notice.

E3194/ET9287/0909

The Haigh Engineering Company Limited

Alton Road, Ross-on-Wye, Herefordshire HR9 5NG

Tel: 01989 763131

Fax: 01989 768777

info@haigh.co.uk



WRAS
APPROVED
PRODUCT



The pulp used for manufacturing this paper is from mills that are certified to ISO 14001. ECF pulp sourced from sustainable forests.