

AMITY OPA RTU

READY TO USE HIGH LEVEL DISINFECTANT / COLD
STERILANT FOR ENDOSCOPES & MEDICAL DEVICES

Amity OPA RTU is safe to use on ALL flexible and rigid endoscopes.



- High Level Disinfectant / Cold Sterilant
- Ready To Use
- No instrument damage
- Excellent material compatibility
- Non-corrosive
- 5 minute soak at 20 °C
- No activation required
- Cost effective
- Tested to EN14885 for disinfectant properties

Concentrate
- 5ltr container

Amity OPA RTU is a superior product for the decontamination of medical, surgical, dental and endoscopic instruments.



DISINFECTANTS

Application

Reprocessing of instruments - high level disinfection

For ALL critical and semi-critical instruments / medical devices including heat sensitive instruments:

For Manual Disinfection of instruments and equipment:


Stage 1:

Make sure that instrument has been thoroughly cleaned (with Viruzyme family or Virudet N) and rinsed.

Stage 2:




1. Fill soaking tray with Amity OPA Ready to use Solution



2. Soak for 5 minutes making sure that all internal channels are flushed through. This ensures High Level Disinfection is attained.

Note: Ensure that all spaces and recesses are completely filled with Amity OPA fluid and no air pockets exist



3. Rinse in sterile water twice and dry with sterile material

The medical device is now ready for use

AMITY OPA

READY TO USE HIGH LEVEL DISINFECTANT /
COLD STERILANT FOR ENDOSCOPES & MEDICAL DEVICES

Amity OPA Microbiological Properties

Active Against	Standards	Contact Time
Spores	EN13704	3 Hours
Mycobacteria	EN14348	5 Minutes
	Mycobacterium terrae	5 Minutes
	Mycobacterium avium	5 Minutes
	EN14563	5 Minutes
Bacteria	Mycobacterium terrae	5 Minutes
	Mycobacterium avium	5 Minutes
Bacteria	EN13727, EN14561	5 Minutes
Viruses	EN 14476 Polio, Adeno and Murine Norovirus	5 Minutes
Fungi / Yeasts	EN14562, EN13624	5 Minutes

Active ingredients: o-Phthalaldehyde (peroxethanoic acid).

CE
1639



Amity Limited, Libra House, West Street,
Worsbrough Dale, Barnsley, South Yorkshire,
S70 5PG, UK

T: +44 (0) 1226 770787

F: +44 (0) 1226 770757

E: sales@amityinternational.com

WWW.AIHC.CO.UK

© Amity Limited all rights reserved 2018

DISTRIBUTED BY: