

Infection Control Equipment for the Hospital Sluice Room



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0000

STANBRIDGE INFECTION CONTROL



Stanbridge Limited have specialised in sluice rooms and infection control equipment since the early 1950s and have pioneered many designs for modern washer/disinfector machines.

Stanbridge Limited have come a long way in promoting infection control, not only in hospitals but also in care homes. Stanbridge are involved with many associations that deal with infection control procedures. We have pioneered the introduction and detail of HTM 20/30, we have had Study Days which have educated those in the care industry to the needs of infection control, and we have sought to bring those in the care sector up to date.

Our comprehensive range of production covers the washer/disinfectors, machines and stainless steel products for sluice rooms, including bespoke work, together with the utensils manufactured by Stanbridge.



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk


Stanbridge reserves the right to make design changes without notification

THE COMPACT SUPERIOR 3

CS3 - 'three utensils at a time' washer/disinfector



Utensil loading configuration



The Compact Superior 3 multi-utensil machine is designed for the full bed occupancy and high dependency nursing homes which have peak loads and a high utensil volume.

The CS3 has the capacity to efficiently wash and disinfect 3 utensils in one operation, cutting water and heating costs and staff waiting time and utensil queueing.

All Stanbridge washer/disinfectors can include the DEOSCALE water treatment system to control limescale build-up and reduce maintenance.

Stanbridge utensil washer disinfectors are manufactured in compliance with the recommendations of British Standard 2745 and to meet the detail requirements of inspectors of the Water Supply Undertakings. Stanbridge Ltd is approved under BS 5750 for quality assurance.

WRAS Approval No. 0310064
Internal break tank has a type 'A' air gap supply



Size	Installation
Height	1095mm (43")
Width	600mm (24")
Depth	600mm (24")
Water Feed	
Supply	Cold
Connection	1/2" BSP male parallel
Consumption	23.5lt (5.25gall) per normal cycle
Overflow	1x19mm (3/4") bore End must be visible (UK requirement)
Waste	82mm (3") 'P' or 'S' trap
Electrical	230V 50Hz 13A single phase Fuse - mains 13A - control 5A
Sound level	73 dB
Washing cycles	1. Normal wash Wash - disinfection at 90°C for 70 seconds - rinse to cool utensils
	2. Super wash Pre flush - wash - disinfection at 90°C for 70 seconds - rinse to cool utensils




Stanbridge Ltd
Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB
Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

THE SLIMLINE

Small size big specification bedpan washer



The "one at a time" traditional bed pan / bottle washer-disinfector is a narrow machine with a small footprint, ideally suited where space in the sluice room is at a premium.

Specifications include choice of hospital or self-generated steam.

A continuous indication, as an option, shows the disinfection temperature. The fail safe interlocks minimise the risk of cross infection, electrical, thermal and mechanical hazards.

Stanbridge now offer Choice - this slim profile machine is one of a range with differing capacity specifications which allow each ward workload to be matched for best infection control and best price, saving space, time and money.

This very high capacity compact machine suitable for multi utensil use suits high bed occupancy and high dependency wards.

Advice to help you with choice is without obligation.

Stanbridge utensil washer disinfectors are manufactured in compliance with the recommendations of British Standard 2745 and to meet the detail requirements of inspectors of the Water Supply Undertakings. Stanbridge Ltd is approved under BS 5750 for quality assurance.

WRAS Approval No. 0310064
Internal break tank has a type 'A' air gap supply



Size

Height
Width
Depth

Installation

1500mm (60")
520mm (21")
490mm (20")

Water Feed

Supply Hot and Cold or Cold
Connection 1/2" BSP male parallel
Minimum Flow 0.19lt/sec (2.65gal/min)
Max Pressure 413kPa (60p.s.i.)
Temperature Cold Ambient
Hot 60-65°C

Overflow

1x19mm (3/4") bore
End must be visible (UK requirement). Pipework by client.

Chamber exhaust

Internal venting system.

Waste

82mm (3") 'P' or 'S' trap

Electrical

230V 50Hz 13A single phase
Fuse - mains 13A - control 5A
Steam Generator Heater 3.0kW
Pump 0.37kW
Consumption 0.21kWh per cycle (average)

Sound level

73 dB

Washing cycles

1. Normal wash cycle
Wash - disinfection at 80°C for 70 seconds - rinse to cool utensils
2. Super wash cycle
Pre flush - wash - disinfection at 80°C for 70 seconds - rinse to cool utensils



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

THE SLIMLINE LOW PROFILE

Low profile washer/disinfector - fits work top units



The Slimline Low Profile multi utensil washer/disinfector has been designed to provide a machine layout with space above where windows and fittings are unavoidable. Worktop areas are preserved. It is ideally suited for surgical wards where the range of utensils required to be disinfected is wide, but daily usage is limited.

Specifications include choice of hospital or self-generated steam.

A continuous indication as an option shows the disinfection temperature. The fail safe interlocks minimise the risk of cross infection, electrical, thermal and mechanical hazards.

Stanbridge now offer Choice - this low profile machine is one of a range with differing capacity specifications which allow each ward workload to be matched for best infection control and best price, saving space, time and money.

Advice to help you with choice is without obligation.

Stanbridge utensil washer disinfectors are manufactured in compliance with the recommendations of British Standard 2745 and to meet the detail requirements of inspectors of the Water Supply Undertakings. Stanbridge Ltd is approved under BS 5750 for quality assurance.

WRAS Approval No. 0310064
Internal break tank has a type 'A' air gap supply



Size	Installation
Height	1000mm (40")
Width	520mm (21")
Depth	600mm (24")

Water Feed	
Supply	Hot and Cold or Cold
Connection	1/2" BSP male parallel
Minimum Flow	0.19lt/sec (2.65gal/min)
Max Pressure	413kPa (60p.s.i.)
Temperature	Cold Ambient
	Hot 60-65°C

Overflow	1x19mm (3/4") bore End must be visible (UK requirement). Pipework by client.
-----------------	---

Chamber Exhaust	Internal venting system.
------------------------	--------------------------

Waste	82mm (3") 'P' or 'S' trap
--------------	---------------------------

Electrical	230V 50Hz 13A single phase Fuse - mains 13A - control 5A Steam Generator Heater 3.0kW Pump 0.37kW Consumption 0.21kWh per cycle (average)
-------------------	---

Sound level	73 dB
--------------------	-------

Washing cycles	1. Normal wash cycle Wash - disinfection at 80°C for 70 seconds - rinse to cool utensils
	2. Super wash cycle Pre flush - wash - disinfection at 80°C for 70 seconds - rinse to cool utensils



Stanbridge Ltd
Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB
Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

DEOSCALE

Water treatment by DeoScale

DeoScale is a water additive to aid cleaning and the build-up of hard scale deposits in automatic machines that wash and disinfect reusable bedpans and commode pots.

Area Scale

In some areas the build-up of limescale is high enough to reduce the effectiveness of the washer/disinfector by causing surface deterioration of utensils. Appearance suffers.

Why water treatment?

Water for washing containers for human waste should be conditioned in the same way as dishwashers for household plates and cups to ensure the fatty residue is removed by reduction in the water surface tension.

The centre of excellence for clinical safety of utensils is hosted by the Department of Health and contained in the recommendations of HTM20/30 for best practice.

Guidelines for the use of washer/disinfectors for handling human waste containers include necessary wash treatments to achieve good outcome.

DeoScale fluid

DeoScale fluid has two components -a non-foaming detergent and a descaler. The addition of a small measured dose of DeoScale fluid to the wash cycle reduces hard water scale build-up. DeoScale is a low-hazard liquid detergent/descaler designed for auto dosing. It is efficient in removing organic and inorganic soils. Low dose in soft water areas reduces water browning.

Background

There are many 'fixed gadget' devices offered to control or treat limescale in water supplies. Confidence in outcome is inconsistent usually because of a dependency on flow rate and limescale levels.

In consequence DeoScale has been developed to control scale build-up in water systems as a new technology for bedpan or commode washer to provide a predictable outcome for varying degrees of district hardness as well as water flow rates. The dose can be set to match local conditions.

Benefits

The system is based upon the accurate metered introduction of the specially formulated fluid into one break tank during each cycle. The chemical composition of this fluid is such that the component of temporary hardness in the water is disorientated and prevented from forming a hard scale layer. Any precipitate passes to drain in each water flow cycle.

New equipment remains pristine. Scaled chambers in old equipment return to 'as new' after a once only hand descale. Existing scale on utensils over a period slowly re-combines and surfaces to 'bright'. Pipe work is cleaned. Organic grubbiness is removed.

Breakdowns and engineers' time has fallen on trial sites as the internal system scale reduces by slowly returning to a soft flushable character.

Cost effectiveness can be demonstrated. Utensil management is excellent and in consequence, so is care staff reaction. Soil removal and control of infection is enhanced.

Based upon the water provider's analysis of the water supplied to the washer/disinfector, the dose required can be calculated or put to work on a 'trial and error' basis.

The count of machine cycles per day for each dirty utility room will identify a budget cost.

Hazard is low - the liquid is classified as a mild irritant.

Equipment

The basic dosing facility can be fitted in the factory for new installations or site added and consists of:

1. Metering pump
2. Fluid carrier for a 5 litre container suitable for wall mounting

DeoScale fluid is packed in 5 litre screw top plastic containers.

Typical sales pack is 4x5 litre boxed suitable for groupage transport.

Contract quantities per year are negotiable. Please contact our Spares Department for prices.

Consumption

The DeoScale consumption for a typical care home or ward with 30 beds would be around 1x5 litres per month for a single machine. In a high scale district this may rise to 2x5 litre per month.



Stanbridge Ltd
Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB
Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

THE GIRVAN RANGE

Easy Clean Utensils



Commode pot

The Stanbridge round commode pot with lid is designed without crevices allowing easy cleaning. With an anti peel surface for longer life.



Commode Insert pan

The insert commode pan improves the commode chair use, as well as sample gathering in the toilet suite.



Slipper pan

The slipper pan has no end pocket to trap soil and eliminates cleaning problems. Three hand grips allow easy positioning. An indent on the plateau reduces point load on the spine. Smooth movement over the mattress is assured.



In Bed pan

The in bed pan gives confidence, stability and secure feeling for sitting patients. The buttress edge gives extra support on the mattress, especially backwards. A raised edge helps location and comfort.



Urine bottle

Robust opaque polythene urine bottle, graduated in ml to 1 litre. Complete with snap lid for transit safety.



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB


Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0006

THE GIRVAN RANGE

Amethyst Racks



The Amethyst Economic (non-heated) quality stainless steel racks are designed to store bedpans, commode pots and urine bottles efficiently and safely. All racks are manufactured to an exceptionally high standard in electro polished 304 stainless steel with a removable stainless drip tray for easy cleaning and disinfection.

The Stanbridge range of washer/disinfectors are an ideal partner for these utensils.

Bedpan rack		order code
Qty of pans	4	Amethyst E BP4
	6	Amethyst E BP6
	8	Amethyst E BP8

Urine Bottle Rack		
Qty of bottles	4	Amethyst E UB4
	6	Amethyst E UB6
	8	Amethyst E UB8

Commode Pot rack		
Qty of pots	4	Amethyst E 4
	6	Amethyst E 6
	8	Amethyst E 8

Urine Bottle and Bedpan rack		
Qty of each	4	Amethyst E UBBP 4
	6	Amethyst E UBBP 6
	8	Amethyst E UBBP 8

Urine Bottle and Commode Pot rack		
Qty of each	4	Amethyst E CP 4
	6	Amethyst E CP 6
	8	Amethyst E CP 8



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0007

THE FORTH RANGE

Stainless steel deep sinks, dedicated handwash bowls, taps and traps



Combined Janitorial and Handwash Unit

Manufactured from grade 304 stainless steel and supplied with adjustable legs, hinged bucket grid, bowl with monobloc swivel mixer tap and combined waste and overflow.

Dimensions: W500 x D750 x H880mm



Free Standing or Wall Mounting Butler Sinks

Manufactured from 1.2mm thick grade 304 satin finish stainless steel with 38mm waste, 2 tapholes, **NO** overflow or plug, 100mm rear upstand and anti-drip edges all round. Mounted on 30mm square tubular stainless steel uprights with adjustable feet but no undershelf. Can be used in conjunction with the 'TAY' range of stainless steel drainers or independently.

Sink 1

Bowl capacity 48 litres W400 x D400 x H300
Overall dimensions: W500 x D600 x H850



Sink 2

Bowl capacity 60 litres W500 x D400 x H300
Overall dimensions: W600 x D600 x H850

Sink 3

Bowl capacity 72 litres W600 x D400 x H300
Overall dimensions: W700 x D600 x H850



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0008

THE FORTH RANGE

Wash Hand Basins



WHB320R

Hand basin constructed of 304 Grade stainless steel. Finish bright/satin. Complete with skirt/apron, bolt set, waste fitting /NO overflow and plug, care leaflet. Taps and traps extra.

Bowl size 320dia x 105mm depth

Overall size: W400 x D350 x H105mm



STHB02

Corner mounted hand basin constructed of 0.9mm thick 304 grade satin finish stainless steel for corner fixing and complete with waste.

Taps and traps extra.

Bowl size: 212 dia x 90mm deep.

Overall size: W305 x D305 x H90mm



BS205

Single wash hand basin manufactured from 0.9mm thick 304 grade satin finish stainless steel with integral soap trays, tap holes and wall brackets.

Bowl size: W500 x D335 x H155mm

Overall size: W550 x D335 x H155mm



STAN5

Wall mounted hand basin manufactured from 1.2mm thick 304 grade bright/satin finish stainless steel. Tap holes 28mm dia. Supplied with skirt, wall fixing.

Note: Taps and traps are extra

Bowl size: W300 x D150 x H100mm

Overall size: W350 x D260 x H185mm



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0009

THE FORTH RANGE

Tap Range and Bottle Traps



STAN EC 04

Mono Bloc mixer tap sink mounted. Chrome plated quarter turn compression valve - 3" lever length (6" lever also available)



STAN ST 04

'H' Section mixer taps sink mounted. Chrome plated quarter turn ceramic - 6" lever length



STAN EC 03

High neck tap sink mounted. Chrome plated quarter turn ceramic valve - 3" lever length. (6" lever also available)



STAN ST 02

Bib taps for wall mounting. Chrome plated quarter turn ceramic valve 3" lever length (6" lever also available)



BOTTLE TRAPS

Chrome plated bottle traps 1 1/4" and 1 1/2"



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

THE TAY RANGE

Ancillary Sluice Room Furniture

Stainless steel work surfaces



Island work surface - free standing

Manufactured from 304 grade stainless steel with 60mm turn down all round and a standard working height of 850mm. Supplied with stand adjustable feet, undershelf and earth tag.

Dimensions: Up to 3500mm long, width 650mm



Sink unit/Work surface

Manufactured from 304 grade stainless steel with 60mm turn down all round and a standard working height of 850mm. Supplied with 'single bowl' and rear upstand.

Dimensions: Up to 3000mm long, width 650mm
Bowl size: 610x457x300mm deep.

Stainless steel drainers



Wall mounted shelf/drainer manufactured from 1.2mm thick grade 304 satin polished stainless steel. Supplied with 2 wall mounting brackets, side aprons, 60mm rear upstand, 2 x 30mm square front support legs with adjustable feet. Butler sink to either lip over left or right.

Type

Drainer 600 - L600mm x W600mm x H960mm

Drainer 900 - L900mm x W600mm x H960mm

Drainer 1200 - L1200mm x W600mm x H960mm

Stainless steel shelves



Heavy duty wall mounted shelves in a standard width of 300mm with integral wall fixing brackets.

Type

Shelving 600 Length 600mm

Shelving 1200 Length 1200mm

Shelving 1500 Length 1500mm

Shelving 1800 Length 1800mm

Shelving 2000 Length 2000mm



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0011

THE THAMES RANGE

Stainless steel sluice units



Garnet Sluice (Cleaner, Sluice/Sink)

Stainless steel cleaner sluice/sink unit.
Stainless steel hinged bucket grid.
Stainless steel wall mounted 9 litre cistern with cover, chrome plated lever action flush handle, side entry supply and overflow and standard cistern kit.
Stainless steel flush pipe.
Upstand to rear and anti-drip edges to front and sides.
Wall mounted wrist lever action hot and cold chrome plated bib taps extended for concealed or exposed supply.
Two cantilever wall brackets.
Hopper 530x530mm.



Sapphire Sluice

Stainless steel cleaner sluice unit.
Stainless steel wall mounted 9 litre cistern with cover, chrome plated lever action flush handle, side entry supply and overflow and standard cistern kit.
Stainless steel flush pipe.
Upstand to rear and anti-drip edges to front and sides.
Surface mounted wrist lever action chrome plated monobloc dual flow mixer tap with swivel spout to sluice.
Two cantilever wall brackets.
Hopper 530x530mm.

THE THAMES RANGE

Stainless steel sluice units



Aquamarine Sluice (Sluice and Work Surface)

Stainless steel left or right handed sluice and work surface.

Stainless steel wall mounted 9 litre cistern with cover, chrome plated lever action flush handle, side entry supply and overflow and standard cistern kit.

Stainless steel flush pipe.

Upstand to rear and anti-drip edges to front and sides.

Surface mounted wrist lever action chrome plated monobloc dual flow mixer tap with swivel spout to sluice.

Two sets of bearers and stainless steel front support legs with adjustable feet or two cantilever wall brackets.

Hopper and work surface 1070x530mm.



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0013

THE THAMES RANGE

Stainless steel sluice units



Opal Sluice (Sluice and Work Surface)

Stainless steel left or right handed sluice and sink unit with work surface. (DUHS)
Stainless steel wall mounted 9 litre cistern with cover, chrome plated lever action flush handle, side entry supply and overflow and standard cistern kit.

Stainless steel flush pipe.

Stainless steel sink to sluice waste connecting pipe.

Upstand to rear and anti-drip edges to front and sides.

Chrome plated 38mm deep seated waste fitting to sink and plastic overflow pipe.

Wrist lever action hot and cold chrome plated monobloc dual flow mixer tap with swivel spout to sluice.

Three sets of bearers and stainless steel front support legs with adjustable feet or three cantilever wall brackets.

Hopper, sink and work surface 1600x530mm.



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0014

THE DERWENT RANGE

Scrub up troughs and sinks

Quartz

(Surgeons' Scrub Up Trough)



Three standard options - 750mm (one station), 1500mm (two station), 2250mm (three station).

Left or right draining with 38mm waste outlet.

Chrome plated 38mm waste fittings.

63mm boxed skirting to allow for the mounting of elbow lever action mixer taps on the trough (taps can also be wall mounted).

Chrome plated bottle tap.

Elbow lever action chrome plated mixer tap with single swivel nozzle for trough (tap holes on 200mm centres) or wall mounting.


Number of taps to suit length of trough.

Full length concealed wall mounting stainless steel joggle bracket.

Earthing terminal.

Super Quartz

(Surgeons' Scrub Up Trough with No-Touch Water and Antiseptic Wash or Soap Dispensing System)



Up to three wash stations each with integral high splash back. 375mm (one station), 750mm (two station), 1500mm (three station), 2250mm.

Left or right draining with 38mm waste outlet.

Chrome plated 38mm waste fittings.

Chrome plated 38mm bottle trap.

Full length concealed wall mounting joggle bracket.

Earthing terminal.

Sensor module (contains sensor, controller, water valve and antiseptic wash/soap pump).

Service module (houses the thermostatic valve and the power pack).

Battery back-up (optional).



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0015

THE DERWENT RANGE

Scrub up troughs and sinks



Zircon Surgeons' Scrub Up Trough

Three standard options - 800mm (one station), 1600mm (two station), 2400mm (three station).

Left or right draining with 38mm waste outlet.

Chrome plated 38mm waste fittings.

63mm boxed skirting to allow for the mounting of elbow lever action mixer taps on the trough (taps can also be wall mounted).

Chrome plated bottle trap.

Elbow lever action chrome plated mixer tap with single swivel nozzle for trough (tap holes on 200mm centres) or wall mounting. Number of taps to suit length of trough.

Full length concealed wall mounting stainless steel joggle bracket.

Earthing terminal.



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0016

THE DERWENT RANGE

Scrub up troughs and sinks



Topaz Plaster Sink and Work Surface

Stainless steel left or right handed plaster sink and work surface.
250mm diameter stainless steel plaster sump and removable mesh filter basket with lifting handle and rebated cover plate.
38mm waste outlet and chrome plated 38mm waste fitting.
Upstand to rear and anti-drip edges to front and sides.
Chrome plated bottle trap.
Elbow lever action chrome plated mixer taps with single swivel nozzle for surface (tap holes on 200mm centres) or wall mounting.
Two sets of bearers and stainless steel front support legs with adjustable feet or two cantilever wall brackets.
Earthing terminal.



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0017

THE TRENT RANGE

Stainless steel cupboards



Cupboard and Sink Unit

Constructed of AISI 304 18/8-10 stainless steel.

Available in 12 models with 2 different basin sizes to cover the vast majority of needs that may arise.

Equipped with slopping groove drainers, so as to guide the draining water into the basin.

Equipped with stainless steel legs which are adjustable in height.



Wall Mounted Cupboard

Constructed of AISI 304 18/8-10 stainless steel.

Available in 4 models of varying sizes.

Equipped with one shelf which is adjustable in height.



Stanbridge Ltd

Leaders in Infection Control Equipment for over 50 Years

Beaumont House, Cray Valley Road, Orpington, Kent BR5 2UB

Tel: 01689 806500 Fax: 01689 806501 email: sales@stanbridge.co.uk www.stanbridge.co.uk

Stanbridge reserves the right to make design changes without notification

H0018

