



UNITED KINGDOM
MEDIVANCE INSTRUMENTS
 Barretts Green Road
 London
 NW10 7AP
 Tel. +44 (0)20 8965 2913
 Fax. +44 (0)20 8963 1270

 www.velopex.com

FRANCE
VELOPEX INTERNATIONAL
 14 route de Montreuil le Gast
 B.P.45
 35520 Melesse
 Tel. +33 (0)2.99.13.22.59
 Fax. +33 (0)2.99.13.22.60

 www.velopex.com

USA
VELOPEX INTERNATIONAL
 105 East 17th Street
 Saint Cloud
 Florida 34769
 Tel. +1 (407) 957 3900
 Fax. +1 (407) 957 3927

 www.velopexusa.com

VISCHECK™

X-RAY QUALITY CONTROL

Don't settle for a Poor X-ray...



- VISCHECK -
Making it
SIMPLE & EASY
for YOU



Yes, please contact me...

For information, prices or order details, just tick the items that interest you, fill in your address and post the card.

- ...for further information.
-I am interested in buying a VELOPEX X-RAY QUALITY CONTROL SYSTEM from you.
- ...to inform me in the future about new products!

Doctor: _____
 Practice: _____
 Address: _____

 City: _____
 Postcode: _____
 Tel. No: _____ Fax. No: _____
 E-Mail: _____
 Preferred Supplier: _____



08.07.08 99 Issue 2



Designed to provide quality assurance and control over x-ray film!

Don't settle for a Poor X-ray - Visccheck it!



Making it **QUICK & EASY** for YOU!


- Tells you when to change your chemicals
- Detects Fog problems on the film
- Helps you in identifying "over temperature" problems
- Helps you in determining whether the film is under or over exposed (e.g. problem with x-ray unit)
- Helps you in determining whether x-ray problems are processor based or whether the problem lies within the x-ray film itself

Use Visccheck to avoid those annoying and time wasting retakes!

This essential aid for x-ray film saves time, money and frustration. With the Visccheck Reader you are able to maintain a record of the quality output from your processor. Helps ensure that x-ray results are correct and accurate every time.



Includes **Electronic Reader** + **20 Visccheck Strips!**



Visccheck enables you to control the conditions under which your x-ray film is being processed and is quick and easy to use!



Simple and easy to understand icons

Chemistry is Excellent
Processor is Excellent




Good enough for processing
- Chemicals will shortly need changing




Chemicals exhausted
- Change chemicals
- Heater temperature too cold



Chemicals are contaminated due to:
- Chemicals too hot
- Too much light in darkroom/daylight loader
- Visccheck film out of date



Simply insert the processed Visccheck Strip into the Visccheck Reader to give you an accurate reading of your processing conditions!



DENTAL Diagnostics