



takara.co.uk/dental

Takara Belmont (UK) Ltd.

London head office & showroom

Belmont House, One St. Andrews Way,
Bow, London E3 3PA UNITED KINGDOM

Telephone +44 (0) 20 7515 0333

Facsimile +44 (0) 20 7987 3596

Email dental@takara.co.uk

Manchester showroom

Telephone +44 (0) 161 743 9992

Showroom visits by appointment



The manufacturer reserves the right to alter specifications without prior notice.



Cleo II
Treatment Centres

Cleo II delivery options :

uniquely ideal :

Cleo II

Equipped for an array of surgery procedures and working methods in complete comfort.

CLEO II 'Surgery System'
Air or Electric version



Advanced technology, improved ergonomics and the hygienic traits of the Cleo II ensure superior control and precision for the dental team.

The flexibility of the Cleo II range ensures that any configuration will equip a dental surgery with the ultimate, complete surgery system. The unique folding and extending chair provides unrivalled comfort and access for both Operator and Patient.

CLEO II 'Flexible System'
chair and cuspidor with Assistant's console



Mobile Cart
Air or Electric version



Cabinet Module
Air or Electric version

- Height-adjustable, extendable, rotating Assistant's console
- Porcelain spittoon
- Rotating, chair-mounted cuspidor (90°)
- Seamless upholstery



Rotating cuspidor unit increases access for an Assistant working closely over the Patient from the 3 o'clock position.

Operator's console is height-adjustable and rotates, permitting precise positioning for ideal working comfort.

'E' version key features ~
LED digital display
Preset micro-motor
Torque control with Auto-Reverse
First-priority handpiece selection

Control and precision are the prerequisites of the Cleo II. The Operator's console rotates smoothly on a compact, base-mount enabling an Operator easy access between the 9 and 1 o'clock positions. Hygiene-conscious at every point, the Cleo II is designed to perform in the most demanding dental surgery environments.

Versatility for **working comfort** :

The Operator can work unobstructed from the front, side or behind the Patient.

Folding legrest and foldaway armrest for perfect Patient entry and exit

Clean and preparation position during Patient consultation improves efficiency



Maximised workflow and a greater variety of operation are the key achievements of the Cleo II treatment centre. An inclined Cleo II chair is a very gentle welcome to any Patient entering the surgery and enables the most pleasant consultation position possible.

PREPARATION

Positioning the Operator's console behind the chair, and instruments out of the Patient's view, reduces a Patient's anxiety and simultaneously offers the perfect clean and preparation arrangement with the Assistant's console.

OPERATION

OPERATION

CONSULTATION

complete, hygienic surgery system :

Cleo II



Surgery System permits multiple dental disciplines with increased Operator comfort.

Ideal for **single-** and **four-handed** working positions.

10

11

12

1

2

3

Extending and rotating the Assistant's console enables very comfortable single-handed operation.

Rotating cuspidor unit increases access for an Assistant working closely over the Patient from the 3 o'clock position.

'E' version

Patient-friendly, synchronised chair :

The unique folding-extending chair provides unrivalled comfort, natural support and enhanced access for both Operator and Patient.



Vacuum solids collector is easily accessed Tubes are easily detached Rotating cuspidor (90°) Foldaway, detachable armrest

Cleo II chair key features ~
 Folding and extending legrest
 Twin-articulating headrest
 Foldaway, detachable armrest
 Synchronised chair movements
 Lumbar support
 Seamless upholstery



Ultra-slim backrest improves Operator's working comfort

Integrated foot control stalks for chair movements including two programmable presets, Last-Position memory and Auto-Return.

The telescopic base is a superior, elegant design solution for easy hygiene maintenance – smooth, white surfaces, with integrated foot controls. Supremely comfortable chair movements, smooth, near-silent and gentle start/stop.

Folding and extending legrest safety feature ensures that the legrest retracts if contact is made with any object.

748mm maximum chair height / 483mm minimum chair height

technically advantaged :

Cleo II

Super ergonomics for Dentist and Patient.

'E' version control panel includes instrument, chair, spittoon and light controls



'A' version control panel includes chair and spittoon controls



Increased range of synchronised chair movements



MOBILE CART (MC) and CABINET MODULE (M) 'A' and 'E' versions



Assistant's console with chair, spittoon and light¹ controls

Control all chair movements by hand or foot: by hand from both the Operator's and Assistant's console; by foot using the stalks built-in to the chair base. Spittoon functions and ON/OFF light control¹ are also activated from both the Operator's and Assistant's console.



820 LIGHT – Unit-mounted or 700 SERIES LIGHT – Unit- or ceiling-mounted

Triple-axis head movement
 White natural beam 5,000 Kelvin
 Touchless sensor ON/OFF control
 Composite Curing Safe Mode

4,200 Kelvin
 Touchless sensor ON/OFF control
 Composite Curing Safe Mode
 Retractable Patient mirror



Operator's foot control for instruments

¹ Light controls on 'E' version only.