



TAKARA BELMONT (UK) LTD.

LONDON HEAD OFFICE
Belmont House, One St. Andrews Way,
Bow, London E3 3PA UNITED KINGDOM
TELEPHONE +44 (0) 20 7515 0333
FACSIMILE +44 (0) 20 7987 3596
EMAIL dental@takara.co.uk
VISIT www.takara.co.uk

MANCHESTER SHOWROOM
TELEPHONE +44 (0) 161 743 9992

SHOWROOM VISITS BY APPOINTMENT



THE MANUFACTURER RESERVES THE RIGHT TO ALTER SPECIFICATIONS WITHOUT PRIOR NOTICE.



VOYAGER II

ADVANCED DENTAL TREATMENT CENTRES

*Voyager II treatment centres are extremely versatile and reliable.
By the nature of its intelligent design and outstanding ergonomics,
it is both a space-saving and efficient working environment.*



The Benchmark For Reliability

LEFT & RIGHT SIDARM DELIVERY

To facilitate absolute versatility for the Dental team.

Voyager II-L enables left and right-handed delivery. The unique, pivoting arm enables the Operator's table to be easily positioned either side of the dental chair.

The changeover can be completed by a single person in just one minute.

ROTATION

Voyager II-L uses a centrally-mounted sidearm to ensure absolute equality for both left and right-handed Operators. The Operator's Table can be positioned where other units simply find it difficult to reach.

ERGONOMICS & HYDRAULICS


The contours and ergonomics of the Voyager Chair are designed to imitate the movements of a Patient's body. Moving parts such as headrests, backrests and bases combine harmoniously in all cases to utilise the best of Belmont's legendary hydraulic technology. The result is a more comfortable and relaxed Patient which in turn aids better Dentistry.



The centrally-pivoting arm makes the Voyager a compact, space-saving unit.



Preparation position



Left-handed



Right-handed

THE VERY VERSATILE VOYAGER

Voyager's high levels of build quality and reliability are supported by free extended warranties.



**OPTIONAL
NDL-VG DENTAL
OPERATING LIGHT**
A smooth, perfectly
balanced and
compact operating
light for ease of use.
A unique sensor
permits touchless
on/off control.
5,000K colour
temperature for
ideal matching
capabilities.




VOYAGER II-L SP FLEXIBLE TREATMENT CENTRE
This unit consists of a rotatable left and right-handed
chair-mounted
cuspidor, available with Mobile Cart or Cabinet Module.



**PERFECTLY BALANCED
SIDEARM DELIVERY**
The sidearm makes
it easy for the
Operator to find
the perfect working
position. The
height-adjustable
table delivers
the instruments
from below, out of
Patient's view.

ASSISTANT'S CONSOLE
The Assistant's Console includes
high vacuum and saliva ejector
outlets. Easily positioned to enable
two and four-handed operation.
A 3-way autoclavable syringe is
available as an optional extra.



HEADREST
The Voyager
chair has a twin
articulating headrest
which allows the
Operator to position
the Patient's head
correctly. It is equally
suitable for both
adults and children.

ERGONOMICS
Synchronised Axis Backrest ensures
optimum Patient comfort at all times. Ultra
slim construction allows ease of access for
both Operator and Assistant. Two pre-set
positions and Last Position Memory make the
Voyager chair simple to operate. A wide range
of movement ensures a very practical working
environment. The minimum seat height is
470mm, and maximum 850mm.

CUSPIDOR
Rotatable (left and
right) Cuspidor unit.
An attractive, hygienic
environment; the
smooth surfaces
and ultra-smooth
cuspidor bowl
allow quick and
efficient, hygienic
cleaning. The (CVS)
solids collector
is conveniently
positioned for easy
cleaning.



UPHOLSTERY
The chair offers super soft upholstery for perfect Patient comfort.
The seamless, easy-clean upholstery offers optimum hygiene levels
and is available in a range of attractive colours.

Voyager II shown in Peppermint seamless upholstery.



VERSATILE, RELIABLE, COMPACT.