



SPECIALISTS IN SUPPLY OF BATTERIES TO THE HEALTHCARE SECTOR

ORIGINAL MANUFACTURER MEDICAL BATTERIES REPLACEMENT MEDICAL BATTERY PACKS
UPS REPLACEMENT BATTERIES EMERGENCY LIGHTING BATTERIES
SEALED LEAD ACID BATTERIES CHARGERS TESTERS



Contents

| | |
|---------|--|
| Page 01 | About Euro Energy Resources Limited |
| Page 02 | Replacement and Original Medical Batteries |
| Page 21 | Alkaline and Rechargeable Batteries/Single Cell Chargers/Lithium Batteries |
| Page 22 | UPS Replacement Batteries |
| Page 23 | Cyclon Sealed Lead Acid Batteries |
| Page 24 | Yuasa Sealed Lead Acid Batteries |
| Page 25 | Panasonic Sealed Lead Acid Batteries |
| Page 26 | Other Sealed Lead Acid Batteries |
| Page 27 | Emergency Lighting Batteries |
| Page 28 | Sealed Lead Acid Battery Chargers/Checkmate III Battery Capacity Tester |
| Page 29 | Battery Terminal and Layout Information |
| Page 30 | Battery Care |
| Page 31 | How to Order and Important Information |





Welcome - About Euro Energy

We are battery specialists and have more than 30 years experience of supplying to the medical sector from our Leicester base. We have established a strong reputation for excellent service, high product quality and a comprehensive product range. This catalogue has been designed to enable you to order a complete range of products from one single source often enabling you to make substantial savings.

Should you not be able to find a battery in our catalogue, have a requirement for an original manufacturer medical battery or you simply wish to discuss any aspect of our product range, then please contact our Customer Service team for further assistance.

There are many benefits in dealing with Euro Energy Resources Limited:

- We specialise in battery supply, concentrating closely on the medical sector - this is an important part of our business
- We continue to expand our range of original manufacturer medical batteries, already covering the majority of the market
- Replacement medical batteries are **CE** marked as legally required by the Medical Devices Directive
- Registered with the Medicines and Healthcare products Regulatory Agency (MHRA) for the manufacture of replacement medical batteries, number CA006924
- BS EN ISO 9001: 2008 Accredited Supplier, ensuring consistent quality
- Where we manufacture, we use identical components to those in the original manufacturer battery, ensuring the same level of reliability and functionality
- Fully compliant with European Medical Devices Directive 93/42/EEC, European Battery Directive 2006/66/EC
- Batteries are tested at various stages of production to ensure quality standards are maintained
- Highest levels of customer care
- Comprehensive stocks ensure many types of battery can be dispatched quickly from stock
- We have acquired a wealth of knowledge and experience through dealing with the medical profession for more than 30 years
- And it's not just medical equipment batteries we supply - see our range of UPS Replacement Batteries (page 22), Emergency Lighting Batteries (page 27) and Chargers & Testers (page 28)
- If a battery you require is not listed, please contact our Customer Service team on (0116) 234 0567 for assistance

Replacement and Original Medical Batteries

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|-------|----------------------|------------|----------|
|-------|----------------------|------------|----------|

3M

| | | | |
|---|-------------|--|----------|
| Surgical Shaver (contact fitting) | 1.2v./1.9Ah | | EE120378 |
| Surgical Shaver (with tags) | 1.2v./1.9Ah | | EE120157 |
| 100, 110, 200, 200A, 210, 270, 285, 300 | 6v./3.8Ah | | 04202 |

ABBOTT

| | | | |
|-------------------------------------|--------------|----|----------|
| TENS Units | 1.2v./150mAh | | 01911 |
| Flexi Flo Portable Nutrient Pump | 6v./1.0Ah | | 04110 |
| Life Care Volumetric Pump 900 | 12v./3.2Ah | | 04129 |
| Flexi Flo Companion | 4v./1.0Ah | | 04207 |
| Life Care Pump 4, 5 (Dryfit) | 8v./2.5Ah | | 04255 |
| Flexiflo-Quantum | 8v./3.2Ah | | 04278 |
| Life Care 5000, 5500, 6000 (Dryfit) | 8v./3.2Ah | | 04278 |
| Flexi Flo Patrol | 6v./4.0Ah | | 04114 |
| Life Care 4200 P.C.A. Pump (Cyclon) | 8v./2.5Ah | CE | EE400130 |
| Life Care 4200 P.C.A. Pump (Dryfit) | 8v./2.5Ah | CE | EE400132 |
| Gemstar Docking Station | 2.4v./3.6Ah | CE | EE090240 |

ACTAMED

| | | | |
|-------------------------|-------------|----|----------|
| 4500 Scout Mini Monitor | 4.8v./1.9Ah | CE | EE120283 |
|-------------------------|-------------|----|----------|

AEROGEN

| | | | |
|---------|--------------|--|----------|
| Aeroneb | 4.8v./450mAh | | EE300230 |
|---------|--------------|--|----------|

AEROSOL PRODUCTS/MG ELECTRIC

| | | | |
|---------------------------|---------------|----|----------|
| SAM 90 Suction Controller | 12v./0.8Ah | CE | EE050030 |
| SAM 418 Suction Unit | 12v./3.0Ah | CE | EE140147 |
| SAM 420 Suction Unit | 2 x 6v./3.0Ah | CE | EE140217 |

AGILENT TECHNOLOGY

(See Hewlett Packard)

AIR SHIELDS (ALSO SEE S & W VICKERS, VICKERS)

| | | | |
|-------------------------|--------------|----|----------|
| 500 Incubator | 7.2v./120mAh | | 03106 |
| Cardiotac Monitor 4000 | 4v./1.0Ah | | 04207 |
| Transilluminator LT1500 | 7.2v./1.9Ah | CE | EE120232 |

ALARIS


(See CareFusion)

ALCON

| | | | |
|-------------------|-------------|----|----------|
| Auto Keratometer | 6v./400mAh | CE | EE040054 |
| Scanner Model 500 | 3.6v./80mAh | CE | EE300010 |

ALOKA


| | | | |
|----------------|--------------|----|----------|
| Scanner SSD650 | 3.6v./300mAh | CE | EE040066 |
|----------------|--------------|----|----------|



REPLACEMENT MEDICAL BATTERIES
3M - BAX


| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|--|------------------------------|--------------------|-----------------|------------------|
| AMEDA | | | | |
| Egnell Elite Breast Pump | 12v./0.8Ah | CE | EE050312 | |
| AND | | | | |
| Ambulatory BP Monitor (TM2421) [insert/rebuild] | 4.8v./0.8Ah | CE | EE050070 | |
| ARCOMED | | | | |
| Volumed VP2000 and 5000 Vol. Infusion Pump | 12v./1.9Ah | CE | EE120122 | |
| Volumed VP2000 and 5000 Vol. Infusion Pump (Ni-MH) | 12v./3.0Ah | CE | EE120366 | |
| VP5005 Micro (Ni-MH) | 12v./2400mAh | CE | EE090165 | |
| VP5005 (Ni-Cd) | 12v./5.0Ah | CE | EE170313 | |
| VP6000, VP7000, SP6000 | 9.6v./2.0Ah | CE | EE080044 | |
| ARJO | | | | |
| Maxilift Hoist | 24v./4.0Ah | CE | 02036 | |
| NDA0100-03 [send for rebuild] | 24v./4.0Ah | | CSD13070 | NEW |
| SARA KPA0100 [send for rebuild] | 24v./5.0Ah | | CSA30060 | NEW |
| ASTROMED | | | | |
| Chart Recorder Dash I, II (S.L.A.) Mark II Recorder | 12v./7.2Ah | | 04274 | |
| ATMOS | | | | |
| Atmovac Suction Pump | 12v./1.2Ah | CE | EE050331 | |
| Atmolit Suction Pump | 12v./1.9Ah | CE | EE120341 | |
| Atmoport Suction Unit | 12v./3.0Ah | CE | EE140289 | |
| B.T. | | | | |
| Claudius II | 7.2v./1.2Ah | CE | EE080018 | |
| BARD | | | | |
| Site Rite 2 [8001C0775] | 14.4v. | | 02135 | Original Battery |
| BV105 | 6v./5.0Ah | | 04021 | |
| 2500 Bladder Scan | 6v./5.0Ah | | 04021 | |
| 2500 Bladder Scan (alternative) | 6v./2.5Ah | | 04011 | |
| SiteRite 3 Ultrasound [insert] | 14.4v./4.2Ah | | EE090206 | |
| BAXTER (ALSO SEE TRAVENOL) | | | | |
| Colleague Infusion Pump [MD12020] | 12v./2.0Ah | | 02124 | Original Battery |
| SMO 200 | 6v./4.0Ah | | 04114 | |
| Flo-Guard 6000, 6100 Infusion Pump | 6v./10.0Ah | | 04118 | |
| 808 Defibrillator | 6v./10.0Ah | | 04118 | |
| Flo-Guard 6201 Infusion Pump | 12v./2.0Ah | | 02124 | |
| Flo-Guard 6300, 6301 Infusion Pump | 12v./3.4Ah | | 04233 | |
| Flo-Guard 6200, 6600 Infusion Pump [2 required per unit] | 4v./6.0Ah | | 04270 | |
| Flo-Guard 100, 200 | 5 x 1.2v./0.8Ah | CE | EE050075 | |
| PAC-X Home Dialysis Unit | 12v./2.5Ah | CE | EE400005 | |
| PAC-T Home Dialysis | 12v./2.4Ah | CE | EE400094 | |
| Colleague Infusion Pump-Memory | 7.2v./2.1Ah | CE | EE460155 | |
| Flo-Guard DCR Pump | 12v./800mAh | CE | EE050176 | |

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|--|----------------------|------------|------------------------|
| BCI INTERNATIONAL | | | |
| Capnocheck Plus | 6v./5.0Ah | CE | EE170444 |
| Capnocheck Plus (SLA) | 12v./2.3Ah | | 04231 |
| 3303 Pulse Oximeter | 9.6v./1.0Ah | CE | EE050389 |
| 4402 Pulse Oximeter, Capnocheck 8200 | 9.6v./1.2Ah | CE | EE050284 |
| BOEHRINGER | | | |
| Coagucheck Plus, System 512 | 7.2v./0.8Ah | CE | EE050306 |
| BRAUN | | | |
| Perfusor FT [flat pack] | 7.2v./1.9Ah | CE | EE120046 |
| Perfusor FT/P [triangular pack] | 7.2v./1.9Ah | CE | EE120093 |
| Perfusor FM Syringe Driver | 7.2v./1.2Ah | CE | EE120311 |
| Infusomat FM Infusion Pump 2000 | 7.2v./2.3Ah | CE | EE140184 |
| Infusomat, Dropmat Secura | 9.6v./4.0Ah | CE | EE170169 |
| FMS Volumetric Pump | 7.2v./1.2Ah | CE | EE120367 |
| BREAS MEDICAL | | | |
| PV401 Ventilator | 12v./6.5Ah | | 04267 |
| PV401, 2 Internal | 24v./1.3Ah | | EE050368 |
| PV401 Memory | 4.8v./0.8Ah | | EE300020 |
| PV403 Internal | 24v./3.6Ah | CE | EE090234 |
| BURDICK/SIEMENS/SPACELABS | | | |
| EK10, Elite E.C.G./II Machine (862278) | 12v./1.9Ah | CE | EE120040 |
| Eclipse 4, 4iEKG, 400, 850, LE11 [send for rebuild] | 16.8v./0.8Ah | | EE050443 |
| CADD | | | |
| Prizm Pump | 12v./2.0Ah | CE | EE080047 |
| CAMBRIDGE HEART INC. | | | |
| VS500 E.C.G. Rec | 2 x 16.8v./1.9Ah | CE | EE120009 |
| CAMBRIDGE INST. (ALSO SEE PICKER INTERNATIONAL) | | | |
| Explorer Mobile X-Ray Unit | 12v./24.0Ah | | 04145 |
| SC1000 E.C.G. Machine [send for rebuild] | 9.6v./0.8Ah | | EE050159 |
| CAMINO LABORATORIES | | | |
| Cranial Press Monitor | 12v./1.2Ah | | 04131 |
| CARDIAC SCIENCE | | | |
| PowerHeart AED G3 PRO [9145-101/9145-001] | 12v./7.5Ah | | 02125 Original Battery |
| PowerHeart AED [9141R-001] | 12v./7.0Ah | | 02134 Original Battery |
| PowerHeart AED G3 [9146-001/9146-102] | 12v./7.5Ah | | 02139 Original Battery |
| PowerHeart AED G3 PRO [9144-001] | | | 02165 Original Battery |
| CARDINAL HEALTH | | | |
| (See CareFusion) | | | |
| CARDIET | | | |
| Model 100 | 12v./0.8Ah | | 04276 |
| AR600 E.C.G. | 9.6v./1.6Ah | | EE050395 |



| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|---|------------------------------|--------------------|-----------------|-----------------------------|
| CARDIO RAPID | | | | |
| Portable E.C.G. Machine | 12v./0.8Ah | | 04276 | |
| CARDIOLINE | | | | |
| Delta 1 and 3 Plus | 24v./0.8Ah | CE | EE050319 | |
| Delta 60 Plus | 24v./1.9Ah | CE | EE120419 | |
| CARDISUNNY | | | | |
| 501A E.C.G. | 9.6v./800mAh | | EE050220 | |
| CAREFUSION | | | | |
| Large Volumetric Pump/GP Pump [1000SP01302] | 9.6v./2.4Ah | | 02126 | Original Battery |
| Asena Syringe Pump GH/CC/PK [1000SP01122] | 7.2v./2.7Ah | | 02133 | Original Battery |
| PC2 Gemini Pump [841026] | 6v./12.0Ah | | 02154 | Original Battery |
| PC1 Gemini Pump [841027] | 6v./7.2Ah | | 02155 | Original Battery |
| Signature Edition Pump [141780] | 12v./1.8Ah | | 02156 | Original Battery |
| 560, 565, 570, 572 Infusion Pump [132412] | 12v./2.5Ah | | 02157 | Original Battery |
| 580, 591, 597, 598 Infusion Pump [129471] | 6v./3.4Ah | | 02158 | Original Battery |
| GW Pump [1000EL00349] | 7.2v./1.25Ah | | 02161 | Original Battery NEW |
| Syringe Driver 741 | 6v./1.2Ah | | 04111 | |
| Syringe Pump P1000, 7000 Series | 6v./3.4Ah | | 04212 | |
| Keo Feed 500 Infusion Pump | 12v./2.1Ah | | 04132 | |
| Keo Feed 3000 Infusion Pump | 12v./2.1Ah | | 04132 | |
| Memory Back-up 770 | 4.8v./50mAh | CE | EE020004 | |
| Thermometer 2000 | 4.8v./0.8Ah | CE | EE050003 | |
| Infusion Pump 630 | 12v./1.9Ah | CE | EE120005 | |
| Inf. Controller 560, 566 [without connector] | 12v./2.5Ah | CE | EE400004 | |
| Inf. Controller 560, 565, 570, 572 [with connector] | 12v./2.5Ah | CE | 400006 | Available from Stock |
| Syringe Pump 700, 710 | 8v./2.5Ah | CE | EE400013 | |
| Syringe Driver 701, 711 | 8v./2.5Ah | CE | EE400042 | |
| Pump 770 | 6v./2.5Ah | CE | EE400065 | |
| Star-Flow Pump 580, 590, 599 Series | 6v./3.4Ah | CE | 400070 | Available from Stock |
| P Series Memory | 2.4v./40mAh | CE | EE300199 | |
| CME | | | | |
| 100P E.C.G. | 12v./0.8Ah | | 04276 | |
| COLIN MEDICAL INSTRUMENTS | | | | |
| Pressmate Blood Pressure Monitor 8800, 8800C, 8800MSP | 12v./2.3Ah | | 04138 | |
| COROMETRICS MEDICAL SYSTEMS | | | | |
| Monitor 555 | 12v./2.3Ah | | 04138 | |
| Monitor 511 | 6v./7.2Ah | | 04273 | |
| COW & GATE | | | | |
| Flocare 500 Enteral Feed Pump | 6v./1.2Ah | | 04111 | |
| Flocare 800 Enteral Feed Pump | 12v./1.2Ah | | 04131 | |

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|--|----------------------|------------|----------|-----------------------------|
| CRITICARE SYSTEM INC. | | | | |
| 504 Pony | 6v./1.2Ah | | 04111 | |
| 503, 1040, 5040 Pulse Oximeter | 6v./1.2Ah | | 04111 | |
| Pulse Oximeter 502, 50245, 504, 50445, 504US, 506, 506DXN | 6v./4.5Ah | | 04294 | |
| 507 Patient Monitor | 12v./2.1Ah | | 04132 | |
| Poet TE602 | 12v./2.1Ah | | 04132 | |
| 504 DX Pulse Oximeter | 6v./3.8Ah | | 04202 | |
| 508 Patient Monitor | 12v./4.0Ah | | 04219 | |
| Poet 600, 601, 602, 602-1, 602-2, 602-3 | 12v./3.4Ah | | 04233 | |
| 500 N.I.B.P. Pulse Oximeter | 12v./1.2Ah | | 04272 | |
| 507NJ NIBP Monitor, 602-14 | 7.2v./2.7Ah | CE | EE090170 | |
| CRITIKON/JOHNSON & JOHNSON/GE HEALTHCARE | | | | |
| Syringe Minder 90 | 6v./1.2Ah | | 04111 | |
| 2100 Volumetric Infusion Pump | 6v./10.0Ah | | 04118 | |
| Dinamap Plus 8700, 8710, 8720, 8725, 9700, 9710, 9720, VSM Plus | 12v./2.3Ah | | 04231 | |
| Compact BP T.S. Monitor | 12v./2.3Ah | | 04231 | |
| Dinamap Pro 100-400 | 12v./2.3Ah | | 04231 | |
| Oxyshuttle Pulse Oximeter [not fast charge] | 9.6v./2.3Ah | CE | EE150020 | |
| Oxyshuttle Pulse Oximeter II | 7.2v./5.0Ah | CE | EE170239 | |
| Dinamap 8100, 8100T, 8110, 8120, 9400, 9300, 9340XL, 9350XL, 9500XL Mon. | 12v./2.8Ah | CE | 400028 | Available from Stock |
| Dinamap ProCare 100, 300, 400 | 6v./3.0Ah | | 04113 | |
| Dinamap Pro 1000 [insert only] | 12v./7.4Ah | | EE090294 | |
| DAMECA | | | | |
| 40099-20 | 14.4v./0.8Ah | | EE050429 | |
| DANBY | | | | |
| (See Baxter) | | | | |
| DATASCOPE/MINDRAY | | | | |
| Accutor+/Trio Monitor [0146-00-0099] | 11.1v./4.4Ah | | 02100 | Original Battery |
| Duo Monitor [0146-00-0079] | 7.4v./6.6Ah | | 02118 | Original Battery |
| Accustat Pulse Oximeter | 12v./2.1Ah | | 04132 | |
| Intra Aorta Balloon Pump 90 [2 required per unit] | 12v./24.0Ah | | 04145 | |
| Intra Aorta Balloon Pump 90 [2 required per unit] | 12v./38.0Ah | | 04146 | |
| Transport Balloon Pump 90L, 90T, 95, 97, 97E [2 required per unit] | 12v./18.0Ah | | 04162 | |
| Passport Transport Monitor | 12v./2.3Ah | | 04231 | |
| Accutor 3, 4 N.I.B.P. Monitor, 3SAT, 4SAT BP | 12v./2.2Ah | CE | EE400072 | |
| DATEX (ALSO SEE OHMEDA MEDICAL PRODUCTS) | | | | |
| TruSat Oximeter [6050-0006-578] | 4.8v./3.2Ah | | 02112 | Original Battery |
| S5 Light Monitor [896895] | 13.2v./1.8Ah | | 02167 | Original Battery NEW |
| AS3-Power Supply | 6v./1.2Ah | | 04111 | |
| Normacap 200 C02 Monitor | 6v./3.0Ah | | 04113 | |
| Satelite Pulse Oximeter [2 required per unit] | 6v./3.0Ah | | 04113 | |
| OSP200 | 6v./3.4Ah | | 04212 | |
| P200 Pulse Oximeter | 12v./3.4Ah | | 04233 | |
| Flexivent II | 12v./2.6Ah | | 04277 | |
| AS5, AS3, LM, LM-1, LM-P1 Compact Monitor | 24v./1.8Ah | CE | EE090107 | |
| AS3 Light Monitor | 13.2v./800mAh | CE | EE050385 | |
| S5 Compact Anaesthesia [2 required per unit] | 12v./3.6Ah | CE | EE090178 | |



| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|---|----------------------|------------|----------|
| AS5, AS3, LM, LM-1, LM-P1 Light Monitor [insert only] | 13.2v./1.9Ah | CE | EE120319 |
| AS5, AS3 Light Monitor [send for rebuild] | 13.2v./1.9Ah | | EE120444 |
| Cardiocap 5 | 12v./2.6Ah | | 04293 |

DATREND SYSTEM

| | | | |
|--------------|-------------|----|----------|
| Oxitest 2000 | 7.2v./1.4Ah | CE | EE090135 |
|--------------|-------------|----|----------|

DEVILBISS HEALTHCARE

| | | | |
|--------------------------------------|------------|----|----------|
| Pulmo Aid Traveller | 12v./3.2Ah | | 04129 |
| Nebulizer | 12v./3.2Ah | | 04129 |
| Compact Suction | 12v./1.3Ah | | 04228 |
| Vac-u-Aid, 7304-A, 7304-D | 12v./1.3Ah | | 04228 |
| Suction Unit 7305 (Sunrise PU101823) | 12v./2.5Ah | CE | EE400159 |

DIAGNOSTIC ULTRASOUND CORP.

| | | | |
|--|-----------|--|-------|
| Bladder Scan BV2500 | 6v./2.5Ah | | 04011 |
| Bladder Scan 500, BVI 2000, BU12000, BZI2500 | 6v./5.0Ah | | 04021 |

DRAEGER

| | | | | | |
|--|--------------|----|----------------|------------------|-----|
| Delta/Gamma Monitor - Internal [5732354] | 14.4v | | 02160 | Original Battery | |
| Oxylog 2000+/3000/3000+ SLA Rechargeable Battery [1843303] | 12v./17.0Ah | | 02162 | Original Battery | NEW |
| Delta/Gamma Monitor - External Lead Acid Battery [5592097] | 12v./2.3Ah | | 02163 | Original Battery | NEW |
| Oxylog 2000+, 3000, 3000+ Ventilator [2M86733] | 11.1v./7.2Ah | | 02164 | Original Battery | |
| Vitolet Monitor | 6v./3.4Ah | | 04212 | | |
| Evita 4 Ventilator (external) | 12v./17.0Ah | | 04235 | | |
| Evita 4 Ventilator (internal) XL, 2 | 12v./3.5Ah | | PLEASE ENQUIRE | | |
| EV800 Ventilator | 12v./2.9Ah | | 04249 | | |
| Julian Ventilator | 3v. | | 50802 | | |
| Evita 4 Ventilator Memory | 3v. | | 50802 | | |
| Oxylog 2000 | 7.2v./1.0Ah | CE | EE050305 | | |
| PEEG Monitor [insert only] | 8.4v./1.9Ah | CE | EE120300 | | |
| 8000 Incubator Alarm | 9v. | | 03104 | | |

DYLADE

| | | | |
|------------------------|------------|----|----------|
| Injectomat 'S' Syringe | 6v./800mAh | CE | EE050122 |
|------------------------|------------|----|----------|

DYNATECH/NEVADA

| | | | |
|-----------------------------------|------------|--|-------|
| Impulse 4000 Defibrillator Tester | 12v./2.1Ah | | 04132 |
|-----------------------------------|------------|--|-------|

EAST RADCLIFFE

| | | | |
|----------------|-------------|--|-------|
| RB2 Respirator | 12v./24.0Ah | | 04145 |
|----------------|-------------|--|-------|

ELMED

| | | | |
|------------------------|------------|--|-------|
| PSO 100 Pulse Oximeter | 12v./2.2Ah | | 04200 |
|------------------------|------------|--|-------|


ENGSTROM

| | | | |
|------------------|----------------------|----|----------|
| Elvira | 1.2/9.6v./270/800mAh | CE | EE050155 |
| Erica Ventilator | 9.6v./1.9Ah | CE | EE120114 |

ESCHMANN

| | | | |
|-----------------|-------------|--|-------|
| Operating Table | 12v./24.0Ah | | 04145 |
|-----------------|-------------|--|-------|

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|--|----------------------|------------|----------|------------------|
| FERNO | | | | |
| Halo Plus Torch | 1.2v./2.2Ah | | 01408 | |
| Streamlight SL20X | 6v./2.3Ah | CE | EE150011 | |
| FIBRE OPTIC MEDICAL | | | | |
| ILS-101 Fibreoptic Illuminator | 7.2v./1.9Ah | CE | EE120372 | |
| FLAEM NUOVA | | | | |
| Mini Transfer Aspirator | 12v./3.2Ah | | 04129 | |
| FRESENIUS | | | | |
| 162 Enteil | 12v./1.3Ah | | 04228 | |
| Frentamat Feeding Pump | 12v./1.3Ah | | 04228 | |
| FUKUDA DENSHI | | | | |
| Cardisuny 501B | 12v./0.8Ah | CE | EE050156 | |
| 5100 | 12v./3.6Ah | CE | EE090103 | |
| Cardisuny 501AXD [no connector] | 12v./1.4Ah | CE | EE090189 | |
| Cardimax FX-2111 E.C.G. Recorder | 9.6v./1.2Ah | CE | EE050398 | |
| GE HEALTHCARE | | | | |
| Dinamap Pro 100-400 BP Monitor [633132] | 12v./2.3Ah | | 02106 | Original Battery |
| Dinamap MPS 722X Monitor [633153/300955] | 12v./3.5Ah | | 02108 | Original Battery |
| Mac 5000/5500 ECG [900770-001] | 18v./3.5Ah | | 02110 | Original Battery |
| Dash 3000/4000 Patient Monitor Battery [2017857-002] | 11.1v./3.52Ah | | 02111 | Original Battery |
| Patient Data Module - PDM [2031069-003] | 11.1v./1.85Ah | | 02149 | Original Battery |
| Vivid I Ultrasound Scanner [2378963-2] | 16.8v. | | 02152 | Original Battery |
| Responder 3000 [30344631] | 12v./2.0Ah | | 02171 | Original Battery |
| MAC 1600 [2035701-001] | 14.4v./2.3Ah | | 02172 | Original Battery |
| MAC VU E.C.G. Monitor [2 required per unit] | 12v./13.0Ah | CE | EE400192 | |
| Tram 7000 Portable Monitor | 12v./2.0Ah | | 04201 | |
| Eagle Monitor-Memory | 3v. | | 50323 | |
| MAC 15 Recorder | 3.6v./0.6Ah | CE | EE070041 | |
| MAC PC, 6, 12 E.C.G. Recorder (3807-001) | 19.2v./1.9Ah | CE | EE120047 | |
| Eagle Monitor-High Capacity | 12v./3.0Ah | CE | EE140099 | |
| MAC 8 Monitor, Max Personal | 19.2v./5.0Ah | CE | EE170242 | |
| Dash Monitor 3000, 4000 | 11.1v. | | 02050 | |
| Cardiosmart Mac 1200 E.C.G. Recorder, Cardiosmart ST | 18v./1.9Ah | CE | EE120345 | |
| 413079-005 [retain backplate] | 13.2v./1.4Ah | CE | EE090224 | |
| Eagle Monitor (Standard Capacity) | 12v./2.5Ah | CE | EE140131 | |
| Eagle 1000 Monitor | 12v./2.5Ah | CE | EE140294 | |
| GOULD | | | | |
| SP1425B Cardiac Output | 6v./1.0Ah | | 04110 | |
| SP2120 Display Monitor [2 required per unit] | 6v./7.0Ah | | 04116 | |
| SP2204 Blood Flow Monitor [2 required per unit] | 6v./7.0Ah | | 04116 | |
| SP1405 Physio Monitor | 12v./7.0Ah | | 04135 | |
| SP2004 Blood Flow Monitor | 12v./7.0Ah | | 04135 | |
| SP2009 Recorder | 12v./7.0Ah | | 04137 | |



| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|---|------------------------------|--------------------|-----------------|-----------------------------|
| GRASEBY | | | | |
| 2000 Series Syringe Pump | 2v./2.5Ah | | 04000 | |
| 3000 Series Syringe Pump [3 required per unit] | 2v./2.5Ah | | 04000 | |
| 915 | 6v./1.2Ah | | 04111 | |
| 3M-AVI200 | 12v./2.1Ah | | 04132 | |
| 304 E.C.G. Monitor | 12v./2.0Ah | | 02124 | |
| 3M-500 Volumetric Pump | 12v./1.3Ah | | 04228 | |
| 9300 AMB PCAS Pump, 903 Oscillomate | 6v./1.4Ah | CE | EE090039 | |
| CADD Prizm Pump | 12v./1.4Ah | | EE090264 | |
| GRAYSTAD | | | | |
| GSI 17 | 7.2v./1.9Ah | CE | EE120324 | |
| GREINER | | | | |
| Labatechnik 240 | 7.2v./120mAh | CE | EE010007 | |
| HAEMOCUE | | | | |
| Haemoglobin Photometer | 6v./0.8Ah | CE | EE050165 | |
| HEMOCHRON (ALSO SEE INTERNATIONAL TECHNADYNE CORP.) | | | | |
| FA 5008 | 7.2v./0.8Ah | | EE050356 | |
| Response | 7.2v./5.0Ah | | EE170460 | |
| 401 Blood Coag. Monitor | 7.2v./1.8Ah | CE | EE120112 | |
| 801 Blood Coag. Monitor | 7.2v./2.3Ah | CE | EE150043 | |
| Junior [send for rebuild] | 8.4v./850mAh | | CSA17050 | |
| HEINE | | | | |
| Ophthalmoscope [X-02.99.382] | 3.5v. | | 02004 | Original Battery |
| Ophthalmoscope [D-01.89.013] | 2.5v./600mAh | | 02089 | Original Battery |
| MPack [X-007.99.676] | 7.2v./4.2Ah | | 02170 | Original Battery NEW |
| Nica Tron N (X-01-99-433) | 2.4v./0.7Ah | | 02022 | |
| Nica Tron S3Z, S2 | 3.6v./0.7Ah | CE | EE130017 | |
| Ophthalmoscope S2Z | 2.4v./0.7Ah | CE | EE130012 | |
| HELLIGE GMBH (ALSO SEE MARQUETTE, HONEYWELL & PPG) | | | | |
| Monitor 300 | 12v./2.3Ah | | 04231 | |
| SMP320 Transport Monitor, SMP300 | 12v./2.3Ah | | 04231 | |
| EK51 | 10.8v./0.8Ah | CE | EE050112 | |
| SCP912 Defibrillator, SCP840 | 18v./0.8Ah | CE | EE050288 | |
| EKG Monitor EK41 | 10.8v./1.9Ah | CE | EE120012 | |
| EK53, EK56, SMS181 | 21.6v./1.9Ah | CE | EE120100 | |
| Serviguard SCP 910, 913, 915, 922 Defibrillator | 12v./1.9Ah | CE | EE120189 | |
| ELD420, ED420, 425, FCP-1 Defibrillator | 12v./3.0Ah | CE | EE140002 | |
| SCP 852 Defibrillator | 12v./2.0Ah | CE | EE140148 | |
| HEWLETT PACKARD/PHILIPS (ALSO SEE PHILIPS) | | | | |
| M3046A Portable Monitor (Veridia M3/M4) ^A | 12v./3.85Ah | | 02025 | |
| 78532 Bedside Monitor | 7.2v./120mAh | | 03106 | |
| 78571, 78572 Recorder | 7.2v./120mAh | | 03106 | |

^A Charger available. Please call our Customer Service Team on (0116) 234 0567 for further information.

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|---|----------------------|------------|----------|
| HEWLETT PACKARD/PHILIPS (ALSO SEE PHILIPS) Continued | | | |
| Central Station | 7.2v./120mAh | | 03106 |
| 1000E Arrhythmia Monitor | 2v./5.0Ah | | 04001 |
| 1504, 1505 E.C.G. Machine | 6v./2.8Ah | | 04112 |
| 78640, Strip Chart Recorder | 6v./2.8Ah | | 04112 |
| M1277A Portable Monitor | 12v./2.3Ah | | 04231 |
| Merlin M1275A Transport Monitor, M1276A | 12v./2.3Ah | | 04231 |
| M1025 Snoopy Monitor | 3v./5.0Ah | | 50902 |
| 78560A Central Station | 3.6v./0.6Ah | CE | EE050237 |
| Codemaster 100, 2475B, M2476B [send for rebuild] | 12v./3.0Ah | | EE140242 |
| 78660, 78670, 78671A Monitor/Defibrillator [insert only] | | | |
| 78350, 78350A, 78353, 78353B Portable | 12v./2.3Ah | CE | EE150003 |
| Neonatal Monitor HP78350, 78353 | 24v./2.3Ah | CE | EE150015 |
| 78619, 78620 Monitor/Defibrillator | 12v./5.0Ah | CE | EE170007 |
| Codemaster 100 M2477B [send for rebuild] | 12v./4.0Ah | | EE170342 |
| 43100, 43110, 43120A, 43130 Monitor/Defibrillator | 12v./2.5Ah | CE | EE400010 |
| Pagewriter M1700, 1702, 1721 E.C.G. Monitor | 18v./7.0Ah | CE | EE400069 |
| Codemaster M1722, M1723, M1724, M1758A, XL, XL-Plus Mon./Defib. | 12v./4.5Ah | CE | 400079 |
| Pagewriter 1770, 200 Series, M2460A, 100 Series | 6v./7.0Ah | CE | 400137 |
| A3 Monitor M2929A | 6v./8.0Ah | | 04031 |
| M1205 | 12v./2.3Ah | | 04231 |
| Opti-Go Ultrasound | 10.8v./7.2Ah | | 02082 |

Available from Stock
Available from Stock

HME

| | | | |
|--------------------------|------------|----|----------|
| Lifepulse 400 | 12v./3.6Ah | CE | EE090243 |
| Lifepulse Mon LP15, LP20 | 12v./2.5Ah | CE | EE140108 |

HONEYWELL (ALSO SEE HELDIGE)

| | | | |
|-------|------------|----|----------|
| CR100 | 12v./1.9Ah | CE | EE120219 |
|-------|------------|----|----------|

HORWELL

| | | | |
|------------------|------------|----|----------|
| Neurothesiometer | 12v./0.8Ah | CE | EE050141 |
|------------------|------------|----|----------|

HOSPAL


| | | | |
|------------------|-------------|----|----------|
| Integra Dialysis | 7.2v./1.9Ah | CE | EE120439 |
|------------------|-------------|----|----------|

HUNTLEIGH NESBITT EVANS

| | | | |
|--|-------------|----|----------|
| Truman Couch Lift | 12v./2.9Ah | | 04249 |
| C880 Contura Bed [4 required per unit] | 12v./1.2Ah | | 04272 |
| Patient Hoist | 12v./2.9Ah | | 04249 |
| 2480 | 9.6v./0.8Ah | CE | EE050098 |
| DFM1 Foetal Monitor | 8.4v./0.8Ah | CE | EE050123 |
| Nimbus II Airbed Memory, Auto Excell | 4.8v./40mAh | CE | EE300158 |
| Evolution Bed | 12v./1.2Ah | | 04131 |
| Baltimore Bed | 6v./1.2Ah | | 04111 |

IMED

| | | | |
|---|-------------|----|----------|
| 800, 900 Series Infusion Pump | 6v./12.0Ah | | 04232 |
| Gemini PC2, PC-2TX Infusion Pump/Controller, Model 1320, 1325 | 6v./12.0Ah | | 04232 |
| Gemini PC4, PC-4-XL (841030) | 12v./3.4Ah | | 04233 |
| Gemini PC1 Infusion Pump/Controller | 6v./7.2Ah | | 04273 |
| Memory Battery PC1, PC2, 3800 Infusion Pumps | 4.8v./80mAh | CE | EE300085 |
| 915 Infusion Pump | 12v./2.5Ah | CE | EE400046 |




| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|---|------------------------------|--------------------|-----------------|-----------------------------|
| IMEX MEDICAL SYSTEMS | | | | |
| Doppler BT0022, 0024 | 12v./300mAh | CE | EE300143 | |
| IMPACT MEDICAL | | | | |
| 306 Suction pump [2 required per unit] | 6v./7.0Ah | | 04116 | |
| INFRASONICS | | | | |
| Adult Star 1010, 1500 | 12v./7.2Ah | | 04274 | |
| Adult Star Ventilator 2000, 1010, 1500 | 12v./7.2Ah | CE | EE400135 | |
| INTERNATIONAL TECHNADYNE CORP. (ALSO SEE HEMOCHRON) | | | | |
| Millenium 35000 | 12v./2.3Ah | | 04231 | |
| Hemachron 401 Blood Coag. Monitor | 7.2v./1.9Ah | CE | EE120112 | |
| Hemachron 801 Blood Coag. Monitor | 7.2v./2.3Ah | CE | EE150043 | |
| Hemachron FA5008 | 7.2v./800mAh | CE | EE050356 | |
| Hemachron Response | 7.2v./5.0Ah | CE | EE170460 | |
| INVACARE | | | | |
| Opale Patient Hoist | 12v./2.9Ah | | 04249 | |
| IN-VIVO RESEARCH INC. | | | | |
| Omega 5600 Pulse Oximeter | 8v./2.5Ah | | 04013 | |
| 4500 Pulse Oximeter | 8v./2.5Ah | | 04013 | |
| Omega 1400, 1445 NIBM | 12v./2.3Ah | | 04138 | |
| Omega 500, 1000, 1100, 1500, 5000 B.P.M. | 12v./7.2Ah | | 04274 | |
| 4500, 4501 Pulse Oximeter (alternative) 4500 Plus, 5600, HB07 | 12v./2.3Ah | | 04138 | |
| 3155A Patient Monitor | 12v./2.3Ah | | 04231 | |
| IVAC/ALARIS | | | | |
| (See CareFusion) | | | | |
| KEELER OPTICAL (SITE MICRO SURGICAL) | | | | |
| Wireless Indirect Ophthalmoscope [1919-P-1013] | 7.4v./1.7Ah | | 02115 | Original Battery |
| Ophthalmoscope [EP39-18918] | 3.75v./2.2Ah | | EE470009 | Original Battery |
| Ophthalmoscope-Vista [length 81mm] | 3.6v./0.7Ah | CE | 100011 | Available from Stock |
| Porta Power | 6v./1.9Ah | CE | EE120110 | |
| Ophthalmoscope 2100 | 2.4v./1.9Ah | CE | EE120149 | |
| Ophthalmic Headlamp | 6v./1.9Ah | CE | EE120330 | |
| Ophthalmoscope-Pre-Vista [length 77mm] | 3.6v./0.7Ah | CE | EE130006 | |
| SmartPack EP39-22079 | 7.2v./2.5Ah | CE | EE090259 | |
| KONTRON (ROCHE) | | | | |
| Vitalogic 4520 [TSGLIIO100298] | 14.8v./5.2Ah | | 02179 | Original Battery NEW |
| 7840, 7845 Pulse Oximeter | 6v./1.0Ah | | 04110 | |
| Cutaneous CO2 Monitor | 6v./10.0Ah | | 04118 | |
| 7650 | 12v./2.1Ah | | 04132 | |
| E.C.G. Monitor 105, 108, 205 | 12v./7.0Ah | | 04137 | |
| 7142 Monitor | 12v./7.0Ah | | 04137 | |
| Microgas 7141, 7142, 7640 Blood Gas Monitor | 12v./7.0Ah | | 04137 | |
| 7141 Monitor | 12v./2.3Ah | | 04138 | |
| 7143 Microrecorder | 12v./2.3Ah | | 04138 | |
| 7341 Micro Recorder | 12v./2.3Ah | | 04138 | |

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|--|----------------------|------------|----------|
| KONTRON (ROCHE) Continued | | | |
| 7501 Defibrillator | 12v./2.3Ah | | 04138 |
| Balloon Pump 7000 [2 required per unit] | 12v./24.0Ah | | 04145 |
| Trakmon 7220 | 12v./5.0Ah | | 04160 |
| Minimon | 12v./2.0Ah | | 04236 |
| LAERDAL | | | |
| Heart Aid 1000 | 12v./1.2Ah | | 04131 |
| Heart Aid 95, Heartstart 1000, 3000 Training Battery | 12v./7.0Ah | | 04137 |
| Heartstart 1000 Defibrillator (training) | 12v./7.0Ah | | 04137 |
| Heartstart 911 | 12v./2.0Ah | | 04201 |
| Heartstart 2000, 3000, 3000QR Defib. (901100) | 12v./2.0Ah | | 04201 |
| Premier Suction Pump | 6v./3.4Ah | | 04212 |
| Suction Pump (Compact, Mini) | 12v./1.3Ah | | 04131 |
| Suction Pump | 2 x 6v./1.9Ah | CE | 120148 |
| Heart Aid 80, 88 Defibrillator | 12v./2.5Ah | CE | EE400109 |
| LSU (780400) | 12v./2.3Ah | | 04231 |
| Pace Aid 53 | 12v./1.3Ah | | 04228 |
| LIBERATOR | | | |
| Alphatalker | 7.2v./1.2Ah | CE | EE080004 |
| Introtalker IT-1 | 7.2v./3.0Ah | CE | EE140191 |
| LIFECARE | | | |
| Freeway Nebulizer | 6v./3.0Ah | | 04113 |
| LIKO | | | |
| Golvo Patient Hoist | 12v./2.9Ah | | 04249 |
| M D E | | | |
| Escort Monitor | 12v./2.3Ah | | 04231 |
| MAQUET | | | |
| Servo I Ventilator [64-87-180] | 12v./4.0Ah | | 02105 |
| Transport Systems | 12v./10.0Ah | | 04248 |
| Beta Theatre Table | 12v./10.0Ah | | 04248 |
| Table Remote Handset | 10.8v./250mAh | CE | EE020040 |
| 1150 | 12v./8.0Ah | | EE400199 |
| MARDIS | | | |
| Orac Communication Aid | 2 x 3.6v./3.0Ah | CE | EE140254 |
| MARQUETTE (ALSO SEE HELLIGE) | | | |
| (See GE Healthcare) | | | |
| MASIMO | | | |
| Radical Hand Held Oximeter [insert] | 4.8v./1.6Ah | | EE050463 |
| MEDELA A.G. | | | |
| Vario | 12v./2.2Ah | | 04200 |

Available from Stock

Original Battery



| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|---|------------------------------|--------------------|-----------------|
| MEDICAL DIAGNOSTIC EQUIPMENT | | | |
| Escort Monitor | 12v./2.3Ah | | 04231 |
| MEDIX INFUSION SYSTEMS | | | |
| Nebulizer Easimist 85301 | 13.2v./0.6Ah | CE | EE070043 |
| Nebulizer M2057 World Trav. | 13.2v./1.9Ah | CE | EE140156 |
| MEDTRONIC (ALSO SEE PHYSIO CONTROL) | | | |
| Cardio Rhythm Atakr | 12v./7.0Ah | | 04137 |
| 540, 550 Blood Pump [2 required per unit] | 12v./7.0Ah | | 04137 |
| MEG | | | |
| Audiometer FF1 | 8.4/16.8v./0.8Ah | CE | EE050099 |
| MERIAN | | | |
| Meri-Cal | 4.8v./0.8Ah | | EE050178 |
| MICRO MEDICAL LTD | | | |
| Microlab 3000 Spirometer | 4.8v./0.8Ah | CE | EE050196 |
| ML 3300, MK3 SER 4000+ | 7.2v./0.6Ah | CE | EE070052 |
| ML 3300, MK4 SER 9000+ | 7.2v./0.6Ah | CE | EE070067 |
| MLA2000, ML 3300, MK2 SER 2000+ | 7.2v./0.6Ah | CE | EE070076 |
| ML 3500, MLA5009, MK5 SER 15000+ | 3.6v./7.2v./0.6Ah | | EE070089 |
| MINOLTA | | | |
| Jaundice Meter 101 | 4.8v./270mAh | CE | EE040003 |
| Pulsox 7 Pulse Oximeter, 7820 Mini Ox. | 7.2v./1.2Ah | CE | EE050051 |
| JM102 | 4.8v./860mAh | CE | EE300166 |
| MOLIFT | | | |
| Hoist | 13.2v./5.0Ah | CE | EE170385 |
| Hoist [insert] | 14.4v./1.5Ah | | EE120380 |
| Quick Raiser [send for rebuild] | 14.4v./1.5Ah | | EE120435 |
| NELCOR PURITAN BENNETT | | | |
| N185 Memory | 3v. | | 03012 |
| Ventilator 7200 Series | 2v./5.0Ah | | 04001 |
| NPB4000 | 6v./8.0Ah | | 04031 |
| 240 Monitor [2 required per unit] | 6v./1.2Ah | | 04111 |
| Pulse Oximeter N1000 [2 batteries required] | 6v./4.0Ah | | 04114 |
| NPB190, 195, 290, 295, 395, 595 | 6v./4.0Ah | | 04114 |
| Monitor N-6000 | 12v./3.4Ah | | 04233 |
| N3000, 3100, 3200 Monitor/Symphony | 12v./2.0Ah | | 02124 |
| Escort 300 Vital Signs Monitor | 12v./2.3Ah | | 04231 |
| Pulse Oximeter N180, N185 | 12v./2.3Ah | | 04231 |
| BP Monitor N3100 | 3v. | | 50402 |
| Pulse Oximeter N-10 Printer | 3v./1.2Ah | | 50900 |
| Pulse Oximeter N10 | 6v./0.8Ah | CE | EE050054 |
| Pulse Oximeter N100 | 12v./2.5Ah | CE | EE400007 |
| Pulse Oximeter N200 (Cyclon) | 6v./2.5Ah | CE | EE400008 |
| Companion 2801 [send wiring loom] | 12v./4.0Ah | CE | EE400088 |
| NPB-70, 75 [send for rebuild] | 7.2v./2.4Ah | | EE090208 |
| N550 Pulse Oximeter | 12v./1.2Ah | | 04121 |

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|--|----------------------|------------|----------|
| NELLCOR PURITAN BENNETT Continued | | | |
| E300 Patient Monitor | 12v./2.3Ah | | 04231 |
| Model 5500 [send for rebuild] | 14.4v./4.0Ah | | EE090256 |
| N200SLA, N1000 [2 batteries required] | 6v./2.4Ah | | 04204 |
| N1000 [2 batteries required] | 6v./4.0Ah | | 04114 |
| N560, N550B | 9.6v./3.6Ah | CE | EE090298 |

| NESBITT & EVANS | | | |
|----------------------------|------------|--|-------|
| Patient Hoist | 6v./1.2Ah | | 04111 |
| Patient Hoist | 12v./2.9Ah | | 04249 |


| NIDEK | | | |
|--------------------------|-----------|----|----------|
| Keratometer KM500 B.P.I. | 6v./1.2Ah | CE | EE080023 |

| NIHON KOHDEN | | | |
|---|------------|----|----------|
| Cardiofax 5101, 5105, 5151 E.C.G. Machine | 12v./2.1Ah | | 04132 |
| Life Scope 6, 6102A, 6105A, 6543, 6543A, 6543L, 6551, 6551A | 12v./2.0Ah | | 04201 |
| Cardiofax 4201 E.C.G. Machine | 12v./2.0Ah | | 04201 |
| Cardiofax 6000 Series E.C.G. Machine | 12v./2.0Ah | | 04201 |
| Cardiolife Tec. 7100A, 7200A, 7300A Monitor | 12v./2.0Ah | | 04201 |
| 8110 Defibrillator/E.C.G. Monitor | 12v./2.0Ah | | 04201 |
| Cardiofax 8830K 3 Channel E.C.G. Machine | 12v./2.3Ah | | 04231 |
| 8023G, 9010, 9020, 9110, 9130, Cardiofax Gem | 12v./2.3Ah | | 04231 |
| Cardiofax 8420K | 12v./2.3Ah | | 04231 |
| Cardiofax 6851K Single Channel | 12v./0.8Ah | CE | EE050275 |
| Cardiolife Mini Defibrillator Tec4150K | 12v./1.1Ah | CE | EE080006 |
| TEC 8250 [refurbish/insert] | 12v./1.9Ah | | EE120255 |
| OLV 1200K | 12v./1.0Ah | CE | EE050418 |

| NONIN | | | |
|---|-------------|----|----------|
| 8600, 8604, 8604D, 8700, 8800, 8804 Series Pulse Oximeter | | | |
| Microspan 1040, 8604F | 6v./1.9Ah | CE | EE120308 |
| 2500 PalmSAT | 4.8v./1.8Ah | CE | EE050375 |
| 9600, 2120 Pulse Oximeter | 7.2v./4.0Ah | | EE090266 |

| NOVAMETRIX MEDICAL | | | |
|--|------------|--|-------|
| 807, 811 Trans Oxygen Monitor | 6v./1.2Ah | | 04111 |
| Nova 500 | 6v./3.0Ah | | 04113 |
| 505, 515, 515A, 520, 520A, 800, 840, 2001 Pulse Oximeter | 12v./2.1Ah | | 04132 |
| 1200, 1260, 7000 CO2 Monitor | 12v./4.0Ah | | 04135 |
| 903 E.C.G. and Apnea Monitor | 12v./4.0Ah | | 04135 |
| 7100, 78100 CO2 Monitor | 12v./2.3Ah | | 04138 |
| Transcutaneous Monitor 840 VFD | 12v./2.3Ah | | 04138 |
| Transcutaneous Monitor 807, 811 | 6v./1.3Ah | | 04226 |

| OHMEDA MEDICAL PRODUCTS (ALSO SEE DATEX) | | | |
|---|--------------|--|-------|
| Infant Warmer 3000, 3300, 3500, 4000 | 7.2v./120mAh | | 03106 |
| Care Plus Incubator | 7.2v./120mAh | | 03106 |
| Syringe Pump 9000 | 6v./1.2Ah | | 04111 |
| Biox. IVA, Oximeter | 12v./2.1Ah | | 04132 |
| Pulse Oximeter 5000 | 12v./7.0Ah | | 04137 |



| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|----------------------------------|------------------------------|--------------------|-----------------|----------------------|
| Pulse Oximeter 3740, 3760, 3760P | 12v./2.2Ah | | 04200 | |
| Modulus 2 Plus Ventilator | 6v./12.0Ah | | 04232 | |
| Monitor 5400, 5410, 5420, 6800 | 4.8v./1.9Ah | CE | EE120053 | |
| Ventilator Model 7800, 7810 | 4.8v./1.9Ah | CE | EE120203 | |
| 3800, 3900 Memory | 4.8v./80mAh | CE | EE300085 | |
| Medipoint 20 | 4.8v./300mAh | CE | EE300117 | |
| Oxycap Monitor 4700 | 12v./2.5Ah | CE | EE400114 | |
| 3700 Pulse Oximeter | 8v./2.5Ah | CE | EE400124 | |
| 3800, 3900, 3901 Pulse Oximeter | 8v./3.2Ah | CE | 400134 | Available from Stock |
| Aestiva 5/7900 | 12v./2.2Ah | CE | EE400145 | |
| Dash Monitor 3000 | 11.1v. | | 02050 | |
| Aestiva 5/7100 | 6v./4.5Ah | | EE400216 | |

OSCAR BOSCAROL

| | | | |
|---------------------|------------|--|-------|
| Suction Pump 032012 | 12v./4.0Ah | | 04135 |
|---------------------|------------|--|-------|

OXFORD

| | | | |
|-----------------|------------|--|-------|
| Hoist | 6v./12.0Ah | | 04120 |
| Hoist Model 125 | 12v./2.9Ah | | 04249 |

PACETECH

| | | | |
|--------------|-----------|--|-------|
| NIPD Monitor | 6v./7.0Ah | | 04116 |
|--------------|-----------|--|-------|

PALCO LABORATORIES

| | | | |
|-------------------------|-------------|----|----------|
| 300, 340 Pulse Oximeter | 7.2v./1.2Ah | CE | EE050326 |
| 400, 500 Pulse Oximeter | 12v./0.8Ah | CE | EE050374 |

PARI

| | | | |
|-----------|------------|----|----------|
| Nebulizer | 12v./1.9Ah | CE | EE120393 |
|-----------|------------|----|----------|

PARKS ELECTRONICS LABS

| | | | |
|---|-------------|----|----------|
| 925L Blood Flow Doppler | 6v./1.0Ah | | 04110 |
| 611, 809, 811L, 811S, 911L, 911S, 915L Blood Flow Doppler | 12v./1.2Ah | | 04131 |
| 1010 Blood Flow Doppler | 12v./7.0Ah | | 04137 |
| 613, 614, 615, 711, 801B, 808, 811, 812, 821, 822, 912 | | | |
| 915 Blood Flow Doppler | 12v./300mAh | CE | EE300003 |

PEGASUS

| | | | |
|------------|------------|--|-------|
| Pro Active | 12v./3.2Ah | | 04129 |
|------------|------------|--|-------|

PENLON BROMSGROVE

| | | | |
|----------------|-------------|----|----------|
| Jet Ventilator | 4.8v./0.8Ah | CE | EE050213 |
|----------------|-------------|----|----------|

PETERS


| | | | |
|----------------------------|-------------------|----|----------|
| Peters Audiometer AP22 | 14.4v./7.2v/0.6Ah | CE | EE070002 |
| Peters Audiometer AP26, 27 | 12v./0.6Ah | CE | EE070003 |

PHILIPS (ALSO SEE HEWLETT PACKARD)

| | | | | |
|--|------------|--|-------|------------------|
| Heartstart XL Defibrillator [M3516A] | 12v./2.0Ah | | 02101 | Original Battery |
| Heartstart FR2+ Defibrillator [M3863A] | 12v./4.2Ah | | 02102 | Original Battery |

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|--|----------------------|------------|----------|-----------------------------|
| PHILIPS (ALSO SEE HEWLETT PACKARD) Continued | | | | |
| IntelliVue MP20/MP30/MP40/MP50 Monitor [M4605A] | 10.8v./6.0Ah | | 02103 | Original Battery |
| Heartstart Forerunner AED [BT1] | 18v./1.3Ah | | 02104 | Original Battery |
| Heartstart MRX Monitor/Defibrillator [M3538A] | 14.8v./6.0Ah | | 02113 | Original Battery |
| Pagewriter Trim I II III EKG [989803130151] | 12v./2.9Ah | | 02116 | Original Battery |
| Intellivue MP2/X2 Monitor [M4607A] | 10.8v./1.0Ah | | 02130 | Original Battery |
| Heartstart FRx/HS1 AED [M5070A] | 9v./4.2Ah | | 02136 | Original Battery |
| Pagewriter Touch EKG [989803129131] | 10.8v./7.2Ah | | 02137 | Original Battery NEW |
| Suresigns VM4, VM6, VM8, VS2, VS3 Monitor [989803144631] | 11.1v./7.2Ah | | 02147 | Original Battery |
| Telemon Monitor M2636A [M4790A/989803109821] | | | 02150 | Original Battery |
| PM9086, 95, 97 Scopemeter | 2 x 2.4v./2.2Ah | CE | EE140172 | |
| A1 Monitor | 6.0v./4.2Ah | | 04282 | |
| A3 Monitor | 6.0v./8.0Ah | | 04031 | |
| PHYSIO CONTROL (ALSO SEE MEDTRONIC) | | | | |
| Lifepak 12 Defibrillator - Sealed Lead Acid [11141-000028] | 12v./2.5Ah | | 02066 | Original Battery |
| Lifepak 1000 Defibrillator [11141-000100] | 12v./6.2Ah | | 02107 | Original Battery |
| Lifepak 500 Monitor/Defibrillator [11141-000147] | 12v./7.5Ah | | 02109 | Original Battery |
| CR+ AED [11403-000002] | 12v./1.5Ah | | 02117 | Original Battery |
| Lifepak 20 Defibrillator/Monitor [11141-000068] | 12v./3.0Ah | | 02120 | Original Battery |
| Lifepak 15 Defibrillator [21330-001176] | 11.1v./5.7Ah | | 02141 | Original Battery NEW |
| Lifepak 20E Monitor/Defibrillator [11141-000112] | 11.1v./6.0Ah | | 02142 | Original Battery |
| Lifepak 12 Defibrillator - Ni-Cd [11141-000149]* | 12v./1.6Ah | | 02159 | Original Battery |
| VSM-3 Cardiac Monitor [2 required per unit] | 6v./1.2Ah | | 04111 | |
| Lifepak 6, 6S and 7 Defibrillator | 12v./1.9Ah | CE | EE120007 | |
| Lifepak 5 Monitor/Defibrillator/Voice recorder [std/fast charge] | 12v./1.0Ah | CE | EE120207 | |
| Lifepak 8 Defibrillator | 2 x 6v./1.9Ah | CE | EE120256 | |
| Lifepak 10 and 11 Monitor/Defibrillator [3 required per unit] | 12v./1.0Ah | CE | EE120305 | |
| Lifepak 6 and 7 Monitor, VSM-2 E.C.G. Monitor | 14.4v./1.9Ah | CE | EE140011 | |
| Lifepak 8 Monitor | 2 x 9.6v./2.3Ah | CE | EE150004 | |
| Lifestat 100 NIBP Monitor | 8v./2.5Ah | CE | EE400017 | |
| Lifepak 9 Monitor/Defibrillator (Cyclon) | 16v./2.5Ah | CE | 400115 | Available from Stock |
| Lifepak 200 Monitor/Defibrillator | 8v./2.5Ah | CE | EE400139 | |
| Lifestat 1600 Pulse Oximeter | 8v./2.5Ah | CE | EE400141 | |
| Lifepak 20 Defibrillator | 12v./3.0Ah | CE | EE120417 | |
| PICKER INTERNATIONAL | | | | |
| (See Cambridge Inst.) | | | | |
| POLYSTAN NORTH AMERICA | | | | |
| PO2 Monitor | 9.6v./5.0Ah | CE | EE170288 | |
| PPG (ALSO SEE HELLIGE) | | | | |
| ELD420, ED 420, 425, FCP-1 Defibrillator | 12v./3.0Ah | CE | EE140002 | |
| PROPAQ/PROTOCOL (ALSO SEE WELCH ALLYN) | | | | |
| 802LTRN Monitor [008-0857-00] | 8.4v./2.0Ah | | 02148 | Original Battery |
| Propaq Memory | 3v. | | 50312 | |
| Propaq 102, 104, 106 Monitor without SP02 | 8v./2.5Ah | CE | EE400098 | |
| Propaq 102, 104, 106, 204, 206 Monitor with SP02, 246CS | 8v./5.0Ah | CE | EE400105 | |

* Replaces Lifepak 12 (Ni-Cd). 2.4Ah Battery. Old part number 02065.



| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|--|------------------------------|--------------------|-----------------|-----------------------------|
| PULMONETICS | | | | |
| LTV1000 Ventilator | 12v./4.5Ah | CE | EE400171 | |
| PURITAN BENNETT | | | | |
| Model 840 Ventilator | 24v./7.0Ah | | EE400214 | |
| QUEST MEDICAL (ALSO SEE TRIMED CORP.) | | | | |
| BCI-AAPLQBC-1108 | 9.6v./1.2Ah | CE | EE050311 | |
| Actester, Coagulation Timer | 10.8v./2.4Ah | CE | EE090128 | |
| RADIOMETER | | | | |
| Oxy-1, 3, 13 Pulse Oximeter | 6v./1.2Ah | | 04111 | |
| TCM3 CO2 Cutaneous Monitor [send for rebuild] | 7.2v./1.9Ah | | EE140118 | |
| REFLOLUX | | | | |
| Model 2 | 6v./0.5Ah | | 50061 | |
| Model 5 | 6v./130mAh | | 50500 | |
| REMEDEI | | | | |
| ENM Pain Relief Device | 8.4v./200mAh | CE | EE300235 | |
| RESPIRONICS | | | | |
| Harmony Ventilator 9303 | 13.2v./10.0Ah | CE | EE180113 | |
| Harmony Battery (charger suitable for 9303 battery [non-original]) | | CE | EE460300 | |
| RIGEL | | | | |
| Multicare Monitor 409-Mempac | 3.6v./150mAh | | 03010 | |
| 304, 309 E.C.G. Monitor | 6v./2.8Ah | | 04112 | |
| S & W VICKERS (ALSO SEE AIR SHIELDS) | | | | |
| Infant Intensive Care C100 Alarm | 7.2v./120mAh | | 03106 | |
| Infant Warmer 786 | 7.2v./120mAh | | 03106 | |
| Diascope Traveller | 12v./2.3Ah | | 04231 | |
| Defi.2 Monitor/Defibrillator [no connector] | 2 x 12v./1.9Ah | CE | EE120020 | |
| DMS 700, 730, 750 Monitor/Defibrillator | 12v./2.3Ah | CE | EE150009 | |
| 2241 | 12v./2.3Ah | CE | EE150044 | |
| Diascope 2 E.C.G. Monitor-2001 | 12v./5.0Ah | CE | EE170034 | |
| C450 Incubator | 8.4v./300mAh | CE | EE300160 | |
| CardioAid LS, MC, MC+, 25, 50, 100, 59439-6 | 12v./1.9Ah | CE | EE120124 | |
| SCALE-TRONIX INC. | | | | |
| Bed Scale 2500 [4 required per unit] | 6v./1.2Ah | | 04111 | |
| Scale 2001, 4004, 4800, 5005, 5700, 6006, 6500, 6700, 20009 | 6v./0.8Ah | CE | EE050007 | |
| SCHILLER | | | | |
| FRED Easy Defibrillator [Easybat] | 12v./2.8Ah | | 02140 | Original Battery |
| FRED Easyport Defibrillator [2.230292] | 12v./750mAh | | 02146 | Original Battery |
| FRED Rechargeable Battery [BUFRNC/U16006] | 12v./1.9Ah | | 02168 | Original Battery NEW |
| AT10+ Lithium [3.920509] | 21.6v. | | 02175 | Original Battery NEW |
| AT1 [4.530024] | 12v. | | 02180 | Original Battery NEW |
| AT2 [4.530023] | 12v./2.6Ah | | 02181 | Original Battery NEW |

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|--|----------------------|------------|----------|
| SCHILLER Continued | | | |
| AT10 EKG Machine [2 required per unit] | 12v./0.8Ah | | 04130 |
| AT-1 | 12v./1.2Ah | | 04131 |
| AT-6, 60 EKG Machine | 12v./2.1Ah | | 04132 |
| AT4 and AT5 EKG Machine | 9.6v./1.2Ah | CE | EE050302 |
| CV3 | 9.6v./1.9Ah | CE | EE120059 |
| SECA | | | |
| Model 727 Scales | 9.6v./0.8Ah | CE | EE050166 |
| Model 944, 954, 644, 964 Scales | 7.2v./0.8Ah | CE | EE050200 |
| Cardiostat 3600 | 12v./1.9Ah | CE | EE120358 |
| Bedscales Model 985 | 7.2v./1.2Ah | CE | EE050388 |
| CT61 E.C.G. | 12v./2.1Ah | | 04132 |
| SC1000 E.C.G. Recorder [send for rebuild] | 9.6v./0.8Ah | | EE050159 |
| CT3000 | 12v./1.2Ah | | 04131 |
| SEWARD MEDICAL | | | |
| 9911 E.C.G. Monitor [2 required per unit] | 6v./2.8Ah | | 04112 |
| AFM210 | 9.6v./0.6Ah | CE | EE050174 |
| 9613 Suction Pump | 12v./3.0Ah | CE | EE140248 |
| SHERWOOD MED./AMERICAN HOSP. SUPPLIES | | | |
| Kangaroo Feed Pumps 200, 220, 300, 330 | 12v./1.9Ah | CE | EE120017 |
| Kangaroo Feed Pumps 324, 224 | 7.2v./1.9Ah | CE | EE120032 |
| Kangaroo Feed Pump 2100 | 3.6v./1.9Ah | CE | EE120125 |
| SEIMED | | | |
| QV2 120 Ultrasound Scanner | 3.6v./80mAh | CE | EE300111 |
| SIEMENS/BURDICK/SPACELABS/IVAC | | | |
| Sirecard 440 E.C.G. Machine | 6v./1.0Ah | | 04110 |
| Cardiostat 31, 31S | 12v./2.2Ah | | 04200 |
| Ventilator 710 | 12v./2.2Ah | | 04200 |
| CardioAid LS, MC, MC+ | 12v./1.9Ah | | EE120124 |
| Cardiostat 3300 | 12v./1.9Ah | CE | EE120228 |
| SC7000, 9000XL Internal | 12v./3.4Ah | | 04233 |
| SC7000, 9000 | 2 x 6v./1.4Ah | | EE090263 |
| Medicare | 6v./4.2Ah | | 04282 |
| Servo SC300 | 12v./2.3Ah | | 04138 |
| SC9000XL, SC7000 External | 12v./2.3Ah | | 04231 |
| SIMED | | | |
| S50, S100, S100E, Oximeter | 8v./2.5Ah | CE | EE400096 |
| SONICAID MEDICAL | | | |
| Foetal Heart Monitor D206, D206F | 9.6v./0.8Ah | CE | EE050011 |
| Foetal Heart Monitor D205 | 3 x 8.4v./300mAh | CE | EE300091 |
| SPACELABS MEDICAL | | | |
| Systems Recorder 2446 | 12v./2.1Ah | | 04132 |
| PCI, PC2 Monitor [2 required per unit] | 12v./2.1Ah | | 04132 |
| 90308 P.C. Express Monitor | 12v./2.3Ah | | 04231 |
| 90309 P.C. Scout Monitor | 12v./2.3Ah | | 04231 |



| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. |
|--|------------------------------|--------------------|-----------------|
| 1st Medic 510, 550, 610 Defibrillator | 12v./2.3Ah | | 04231 |
| 90305/6, 90315/6 PC2 Monitor [3 required per unit] | 8v./3.5Ah | | 04242 |
| Power Brick | 12v./0.8Ah | | 04276 |
| 90303, 90303B, 90311, 90312, 90312A, 90315 Memory | 3v./170mAh | | 50800 |
| 90600 Series Clock | 3v./170mAh | | 50800 |
| 90207 | 3v./170mAh | | 50800 |
| Alpha PC1 Monitor (90303/90311/90312B) | 12v./2.5Ah | CE | EE400123 |
| Ultraview 146-0053-00 | 12v./0.7Ah | CE | EE020033 |
| 1050 (Ni-MH) 040-1073-00 | 12v./2.4Ah | CE | EE090241 |
| 1050 (SLA) | 12v./2.3Ah | | 04231 |
| Twin Charger for 1050 Ni-MH Battery | | | 04929 |

SUZUKEN KENZ

| | | | |
|--------------------------|--------------|----|----------|
| 103, 106 E.C.G. Recorder | 21.6v./0.8Ah | CE | EE050279 |
|--------------------------|--------------|----|----------|

T.O.P.

| | | | |
|----------------------|-------------|----|----------|
| Syringe Pump 5100 RD | 9.6v./1.1Ah | CE | EE090035 |
| Syringe Pump 3100 | 12v./1.4Ah | CE | EE090134 |

TERUMO

| | | | |
|----------------------|-------------|----|----------|
| TE-172 Infusion Pump | 9.6v./1.2Ah | CE | EE120386 |
|----------------------|-------------|----|----------|

TOBY CHURCHILL

| | | | |
|----------------------|--------------|----|----------|
| Light Writer SL1 | 1.2v./150mAh | | 03036 |
| CPM Elbow Unit SC104 | 7.2v./3.0Ah | CE | EE140281 |
| Lightwriter | 4.8v./1.4Ah | CE | EE090052 |
| Light Writer SL4 | 4.8v./5.0Ah | CE | EE170391 |
| Lightwriter SL30, 35 | 4.8v./2.1Ah | CE | EE050372 |
| Lightwriter SL86 | 4.8v./5.0Ah | CE | EE170450 |

TORONTO

(See Toby Churchill)

TRAVENOL

(See Baxter)

TRIMED CORP. (ALSO SEE QUEST MEDICAL)

| | | | |
|--|--------------|----|----------|
| Actester, Blood Coagulation Tester PA-1000 | 10.8v./1.9Ah | CE | EE140146 |
|--|--------------|----|----------|

VERATHON MEDICAL

| | | | | |
|---------------------------------|-------------|--|-------|------------------|
| BVI3000 Bladderscan [0400-0039] | 7.2v./1.2Ah | | 02131 | Original Battery |
|---------------------------------|-------------|--|-------|------------------|

VICKERS (ALSO SEE AIR SHIELDS)

| | | | |
|----------------------------------|---------------|----|----------|
| IP5 Syringe Pump | 12v./2.1Ah | | 04132 |
| 77 Transport Incubator | 12v./38.0Ah | | 04214 |
| Transport Ventilator 77, 90 | 12v./0.8Ah | CE | EE050020 |
| Trailblazer Incubator/Ventilator | 7.2v./1400mAh | CE | EE090069 |

VIOMEDIX


| | | | |
|-----------------|-----------|--|-------|
| Feed Pump X700F | 6v./1.3Ah | | 04226 |
| VX600D | 6v./1.2Ah | | 04111 |

VITALOG MONITORING

| | | | |
|--------------------|-------------|----|----------|
| 5000 Sleep Monitor | 7.2v./1.4Ah | CE | EE070042 |
|--------------------|-------------|----|----------|

| MODEL | VOLTAGE/ CAPACITY | CE MARK | PART NO. | |
|---|----------------------|------------|----------|-----------------------------|
| VITALOGRAPH | | | | |
| Spirometer 2120 | 3.6v./0.8Ah | | EE050281 | |
| WEINMANN | | | | |
| Accuvac Rescue | 12v./3.0Ah | CE | EE140288 | |
| WELCH ALLYN (ALSO SEE PROPAQ/PROTOCOL) | | | | |
| AED10 [00185-2] | 12v./2.6Ah | | 02138 | Original Battery |
| AED20 - Non-Rechargeable [00183-OE] | 12v./5.2Ah | | 02143 | Original Battery |
| LXi Spot Vital Signs Monitor [4500-84] | 6v./6.0Ah | | 02145 | Original Battery |
| CP100/CP200 [CP200] | 6v./5.0Ah | | 02166 | Original Battery NEW |
| 62000 Atlas | 6v./7.2Ah | | 04273 | |
| 72900 Microtymp 2 | 2.4v./0.6Ah | CE | EE070054 | |
| 72240 Indirect Ophthalmoscope | 4.8v./1.9Ah | CE | EE120244 | |
| 71670, 72200, 72207 Halogen Ophthalmoscope | 3.6v./0.7Ah | CE | 130005 | Available from Stock |
| 72300, 71000-A, 71020-A, 71000C, 71020-C, 71055-C | 3.6v./0.7Ah | CE | EE130007 | |
| 72000 Pocket Scope | 2.4v./0.7Ah | CE | EE130011 | |
| 5200 Vital Signs Monitor | 6v./4.5Ah | CE | 400161 | Available from Stock |
| 53000 Vital Signs Monitor, 300 Series, 5300 | 6v./4.0Ah | CE | EE400226 | |
| 72250 | 3.6v./1.9Ah | CE | EE120391 | |
| 420, 4200B, 420TB, 42MOB, 42NOB, 42MTB, 42NTB Series | 6v./4.5Ah | CE | EE400220 | |
| 72100, 72110 | 2.4v./700mAh | CE | EE130008 | |
| CP100/CP200 | 6v./5.0Ah | CE | EE400233 | |
| WELMED | | | | |
| P1000, 2000, 3000, 4000 Syringe Pump | 6v./2.8Ah | | 04112 | |
| ZOLL | | | | |
| PD1400/1600/1700/2000/4410/M | | | | |
| Series AED Pro Defibrillator [8000-0299-01] | 10v./2.5Ah | | 02119 | Original Battery |
| XL Smart Defibrillator - With Fuel Gauge [8000-0687-01] | 10v./4.5Ah | | 02128 | Original Battery |
| XL Smart Ready Defibrillator - No Fuel Gauge [8000-0500-01] | 10v./4.5Ah | | 02129 | Original Battery |
| AED PRO R Series Surepower [8019-0535-01] | 10.8v./5.8Ah | | 02153 | Original Battery |
| Auto Pulse [8700-0702-01] | 32.4v./3.2Ah | | 02169 | Original Battery NEW |
| PD1400, 1600, 2000, M series Monitor/Defibrillator | 10v./2.5Ah | | 02028 | |

**IF YOU REQUIRE A BATTERY THAT IS NOT LISTED,
PLEASE CALL OUR CUSTOMER SERVICE TEAM ON
(0116) 234 0567**



Alkaline Batteries **Panasonic**

| MODEL | SIZE | DIAMETER (mm) | HEIGHT (mm) | VOLTAGE | PART NO. | MODEL | SIZE | DIAMETER (mm) | HEIGHT (mm) | VOLTAGE | PART NO. |
|-------|------|---------------|-------------|---------|----------|-------|------|---------------|-------------|---------|----------|
| LR03 | AAA | 10.5 | 44.5 | 1.5v. | 50003 | LR20 | D | 34.2 | 61.5 | 1.5v. | 50000 |
| LR6 | AA | 14.5 | 50.5 | 1.5v. | 50002 | 6LR61 | PP3 | - | - | 9.0v. | 50100 |
| LR14 | C | 26.2 | 50.5 | 1.5v. | 50001 | | | | | | |

PLEASE NOTE THAT PACK SIZES MAY VARY FROM TIME TO TIME. PLEASE CHECK UPON ORDERING. PANASONIC BATTERIES ARE NORMALLY SUPPLIED, BUT BRANDS SUPPLIED MAY CHANGE FROM TIME TO TIME ACCORDING TO AVAILABILITY.

Rechargeable Batteries

| MODEL | CHEMISTRY | SIZE | DIAMETER (mm) | HEIGHT (mm) | VOLTAGE/CAPACITY | PART NO. |
|-------------|-----------|------|---------------|-------------|------------------|----------|
| GP100 | Ni-MH | AAA | 10.5 | 44.5 | 1.2v./0.97Ah | 00218 |
| GP210 | Ni-MH | AA | 26.0 | 49.0 | 1.2v./2.10Ah | 00539 |
| GP350 | Ni-MH | C | 26.0 | 49.0 | 1.2v./3.50Ah | 01423 |
| GP700 | Ni-MH | D | 33.0 | 61.0 | 1.2v./7.00Ah | 01702 |
| GP17R8H | Ni-MH | PP3 | - | - | 8.4v. | 03114 |
| P18N | Ni-Cd | N | - | 10.0 | 1.2v. | 01911 |
| VL2320/IVC* | - | Coin | 23.0 | - | 3.0v. | 03012 |

*with pins

Single Cell Chargers for Rechargeable Batteries

Part No: 04912

- Includes 4 x Ansmann 2400 Batteries
- Suitable for Ni-MH
- Safety timer
- Fast charges 1-4 AA/AAA
- Plug-in type
- Charges PP3 batteries (70mA charge)
- Detects when cells are fully charged

Lithium Batteries

| TYPE | VOLTAGE | DIAMETER (mm) | HEIGHT (mm) | PART NO. | TYPE | VOLTAGE | DIAMETER (mm) | HEIGHT (mm) | PART NO. |
|-------------------|---------|---------------|-------------|----------|--------|---------|---------------|-------------|----------|
| CR1025 | 3v. | 10.0 | 2.5 | 50301 | CR2032 | 3v. | 20.0 | 3.2 | 50312 |
| CR1216 | 3v. | 12.5 | 1.6 | 50303 | CR2320 | 3v. | 23.0 | 2.0 | 50313 |
| CR1220 | 3v. | 12.5 | 2.0 | 50304 | CR2430 | 3v. | 24.0 | 4.0 | 50317 |
| CR1616 | 3v. | 16.0 | 1.6 | 50306 | CR2450 | 3v. | 25.0 | 5.0 | 50318 |
| CR1632 | 3v. | 16.0 | 3.2 | 50308 | CR2330 | 3v. | 23.0 | 3.0 | 50319 |
| CR2016 | 3v. | 20.0 | 1.6 | 50310 | CR2477 | 3v. | 24.0 | 7.7 | 50321 |
| CR2025 | 3v. | 20.0 | 2.5 | 50311 | CR3032 | 3v. | 30.0 | 3.2 | 50322 |
| Ultralife U9VL FP | 9v. | | | 50110 | | | | | |

MANY OTHER BATTERIES ARE AVAILABLE FOR DIFFERENT APPLICATIONS - LIGHTING AND LAMPS - HEARING AIDS - PHOTOGRAPHIC - MEMORY BACK-UP.

PLEASE CALL OUR CUSTOMER SERVICE TEAM ON (0116) 234 0567 FOR FURTHER INFORMATION ON OUR COMPLETE RANGE OF BATTERIES.

UPS Replacement Batteries

| MODEL | PART NO. | MODEL | PART NO. |
|-------|----------|-------|----------|
| RBC1 | EE400257 | RBC24 | EE400270 |
| RBC2 | EE400258 | RBC25 | EE400283 |
| RBC3 | EE400259 | RBC26 | EE400276 |
| RBC4 | EE400260 | RBC27 | EE400277 |
| RBC5 | EE400250 | RBC29 | EE400278 |
| RBC6 | EE400251 | RBC30 | EE400279 |
| RBC7 | EE400252 | RBC31 | EE400270 |
| RBC8 | EE400253 | RBC32 | EE400265 |
| RBC9 | EE400254 | RBC33 | EE400281 |
| RBC10 | EE400261 | RBC35 | EE400282 |
| RBC11 | EE400255 | RBC40 | EE400286 |
| RBC12 | EE400256 | RBC43 | EE400271 |
| RBC13 | EE400268 | RBC44 | EE400287 |
| RBC14 | EE400268 | RBC47 | EE400289 |
| RBC17 | EE400274 | RBC48 | EE400266 |
| RBC18 | EE400269 | RBC54 | EE400270 |
| RBC21 | EE400275 | RBC55 | EE400255 |
| RBC22 | EE400262 | RBC59 | EE400270 |
| RBC23 | EE400264 | RBC62 | EE400263 |



RBC9 EE400254



RBC8 EE400253



RBC12 EE400256



RBC26 EE400276


We would like to point out that Euro Energy is not an authorised APC dealer and this product is not approved by APC.

APC & RBC are registered trademarks of American Power Conversion Corporation. Euro Energy use the above for reference purposes only.

PLEASE NOTE THAT MANY OTHER UPS REPLACEMENT BATTERIES ARE AVAILABLE. PLEASE CALL OUR CUSTOMER SERVICE TEAM ON (0116) 234 0567 FOR FURTHER INFORMATION.

RBC1 EE400257





Sealed Lead Acid Batteries



| CYCLON CELL TYPE | MANUFACTURER PART NO. | DIMENSIONS (mm) | HEIGHT (mm) | HEIGHT TO TERMINALS (mm) | VOLTAGE/CAPACITY | WEIGHT APPROX. (kg) | PART NO. | TERMINAL SIZE (mm) |
|------------------|-----------------------|-----------------|-------------|--------------------------|------------------|---------------------|----------|---------------------|
| D Cell | 0810-0004 | 34 (dia) | 61.0 | 68.0 | 2v./2.5Ah | 0.20 | 04000 | 4.75 x 0.64 Tabs |
| DT Cell | 0860-0004 | 34 (dia) | 96.0 | 103.0 | 2v./4.5Ah | 0.30 | 04006 | 4.35 x 0.64 Tabs |
| X Cell | 0800-0004 | 44 (dia) | 73.0 | 82.0 | 2v./5.0Ah | 0.40 | 04001 | 6.35 x 0.64 Tabs |
| E Cell | 0850-0004 | 44 (dia) | 100.0 | 109.0 | 2v./8.0Ah | 0.60 | 04005 | 6.35 x 0.64 Tabs |
| BC Cell | 0820-0004 | 65 (dia) | 169.0 | 174.0 | 2v./25.0Ah | 1.70 | 04003 | M8+M6 Threaded Post |
| D Monobloc | 0819-0010 | 80 x 46 (LxW) | 70.0 | N/A | 4v./2.5Ah | 0.40 | 04010 | 4.75 x 0.64 Tabs |
| X Monobloc | 0809-0010 | 97 x 54 (LxW) | 77.0 | N/A | 4v./5.0Ah | 0.80 | 04020 | 6.35 x 0.64 Tabs |
| E Monobloc | 0859-0010 | 97 x 54 (LxW) | 102.0 | N/A | 4v./8.0Ah | 1.00 | 04036 | 6.35 x 0.64 Tabs |
| D Monobloc | 0819-0012 | 114 x 46 (LxW) | 70.0 | N/A | 6v./2.5Ah | 0.50 | 04011 | 4.75 x 0.64 Tabs |
| X Monobloc | 0809-0012 | 139 x 54 (LxW) | 77.0 | N/A | 6v./5.0Ah | 1.00 | 04021 | 6.35 x 0.64 Tabs |
| E Monobloc | 0859-0012 | 139 x 54 (LxW) | 102.0 | N/A | 6v./8.0Ah | 1.50 | 04031 | 6.35 x 0.64 Tabs |

Features

- Charging:** To 95% capacity within 1 hour. No requirement to limit charge current.
- Storage Life:** 2 years from fully charged. Check voltage after 6 months.
- Cycle Life:** Up to 2500 cycles @ 30% depth of discharge.
- Float Life:** 15 years at 20°C.
- High current performance:** Up to 26c discharge rate.
- Strength:** Very high resistance to shock and vibration.
- Installation:** Can be used in any orientation.
- Safety:** No loose electrolyte - will not leak acid.
- Transportation:** Non-hazardous, non-spillable cargo classification by land, air or sea.
- Operating temperature:** -65 to +80°C.

PLEASE NOTE THAT OTHER MONOBLOCKS ARE AVAILABLE. PLEASE CONTACT OUR CUSTOMER SERVICE TEAM ON (0116) 234 0567 FOR FURTHER INFORMATION.

Sealed Lead Acid Batteries



NP38-12
PART No: 04141



NP7-12
PART No: 04137



NP4-6
PART No: 04114



NP2.3-12
PART No: 04138

| MODEL | L x W x H (mm) | HEIGHT TO TERMINALS (mm) | VOLTAGE/ CAPACITY | WEIGHT APPROX. (kg) | PART NO. | LAYOUT | TERMINAL |
|-----------------------|-----------------------|--------------------------|-------------------|---------------------|-----------------------|--------|----------|
| NP1-6 | 51.0 x 42.5 x 51.0 | 54.5 | 6v./1.0Ah | 0.25 | 04110 | 5 | A |
| NP1.2-6 | 97.0 x 25.0 x 51.0 | 54.5 | 6v./1.2Ah | 0.30 | 04111 | 1 | A |
| NP2.8-6 | 134.0 x 34.0 x 60.5 | 64.0 | 6v./2.8Ah | 0.66 | 04112 | 1 | A |
| NP4-6 | 70.0 x 47.0 x 102.0 | 105.5 | 6v./4.0Ah | 0.87 | 04114 | 5 | A |
| NP7-6 | 151.0 x 34.0 x 94.0 | 97.5 | 6v./7.0Ah | 1.32 | 04116 | 1 | A |
| NP10-6 | 151.0 x 50.0 x 94.0 | 97.5 | 6v./10.0Ah | 1.99 | 04118 | 1 | A |
| NP12-6 | 151.0 x 50.0 x 93.5 | 97.5 | 6v./12.0Ah | 2.05 | 04120 | 1 | D |
| NP0.8-12 | 96.0 x 25.0 x 61.5 | N/A | 12v./0.8Ah | 0.37 | 04130 | 6 | I |
| NP1.2-12 | 97.0 x 48.0 x 51.0 | 54.5 | 12v./1.2Ah | 0.64 | 04131 | 3 | A |
| NP2.1-12 | 178.0 x 34.0 x 60.5 | 64.0 | 12v./2.1Ah | 0.91 | 04132 | 1 | A |
| NP2.3-12 | 178.0 x 34.0 x 60.5 | 64.0 | 12v./2.3Ah | 0.95 | 04138 | 1 | A |
| NP2.8-12 | 134.0 x 67.0 x 60.5 | 64.0 | 12v./2.8Ah | 1.32 | 04134 | 3 | A |
| NP2.9-12 ^A | 79.0 x 55.5 x 99.0 | 102.0 | 12v./2.9Ah | 1.17 | Now 04249 (PB02.9-12) | 2 | A |
| NP3.2-12 | 134.0 x 67.0 x 60.5 | 64.0 | 12v./3.2Ah | 1.31 | 04129 | 3 | A |
| NP4-12 | 90.0 x 70.0 x 102.5 | 106.0 | 12v./4.0Ah | 1.64 | 04135 | 1 | A |
| NPH5-12 | 90.0 x 70.0 x 102.0 | 106.0 | 12v./5.0Ah | 2.00 | 04160 | 1 | D |
| NP7-12 | 151.0 x 65.0 x 94.0 | 97.5 | 12v./7.0Ah | 2.65 | 04137 | 4 | A |
| NP12-12 | 151.0 x 98.0 x 93.5 | 97.5 | 12v./12.0Ah | 4.00 | 04128 | 4 | D |
| NP17-12 | 181.0 x 76.0 x 167.0 | N/A | 12v./17.0Ah | 6.06 | 04139 | 2 | E |
| NP18-12B | 181.0 x 76.0 x 167.0 | N/A | 12v./18.0Ah | 6.00 | 04162 | 2 | E |
| NP24-12 | 175.0 x 166.0 x 125.0 | N/A | 12v./24.0Ah | 9.68 | 04140 | 2 | J*/E |
| NP38-12 | 197.0 x 165.0 x 170.0 | N/A | 12v./38.0Ah | 13.88 | 04141 | 2 | J*/F |
| NP65-12 | 350.0 x 166.0 x 174.0 | N/A | 12v./65.0Ah | 22.82 | 04142 | 2 | K/G |
| NPC17-12 | 181.0 x 76.0 x 167.0 | N/A | 12v./17.0Ah | 6.36 | 04144 | 2 | E |
| NPC24-12 | 175.0 x 166.0 x 125.0 | N/A | 12v./24.0Ah | 9.64 | 04145 | 2 | J*/E |
| NPC38-12 | 197.0 x 165.0 x 170.0 | N/A | 12v./38.0Ah | 13.95 | 04146 | 2 | F |
| NPC65-12 | 350.0 x 166.0 x 174.0 | N/A | 12v./65.0Ah | 23.00 | 04153 | 2 | K |
| YC33-12 | 195.0 x 129.0 x 179.0 | | 12v./33.0Ah | 11.95 | 04150 | | |
| YC50-12 | 228.0 x 138.0 x 238.0 | | 12v./50.0Ah | 18.71 | 04151 | | |
| YC75-12 | 258.0 x 168.0 x 228.0 | | 12v./75.0Ah | 25.20 | 04152 | | |
| YC100-12 | 344.0 x 171.0 x 280.0 | | 12v./100.0Ah | 39.00 | 04157 | | |

^A Obsolete - an alternative battery (non Yuasa) is available. See page 26. Use part numbers to cross reference details.

* Available with 'L' shaped bracket.

YUASA NP RANGE IS RECOMMENDED FOR GENERAL USE, NPC AND YC RANGES ARE RECOMMENDED FOR CYCLIC USE. PLEASE ENSURE YOU CAREFULLY CHECK BATTERY SPECIFICATIONS AND THE LAYOUT AND TERMINAL INFORMATION ON PAGE 29. ALL BATTERIES HIGHLIGHTED ABOVE ARE OF THE ABSORBENT GLASS MATT (AGM) TYPE. **FURTHER YUASA BATTERIES ARE AVAILABLE, PLEASE CONTACT OUR CUSTOMER SERVICE TEAM ON (0116) 234 0567 FOR FURTHER INFORMATION.**

Sealed Lead Acid Batteries

Panasonic®



LCR063R4 PART No: 04212



LCR064R5 PART No: 04282



LCR121R3PG PART No: 04228



LCR127R2PG PART No: 04274



LCRA1212P PART No: 04211

| MODEL | L x W x H (mm) | HEIGHT TO TERMINALS (mm) | VOLTAGE/CAPACITY | WEIGHT APPROX. (kg) | PART NO. | LAYOUT | TERMINAL |
|---------------------------|-----------------------|--------------------------|------------------|---------------------|------------------------|--------|-------------------|
| LCR061R3PG | 97.0 x 24.0 x 50.0 | 55.0 | 6v./1.3Ah | 0.30 | 04226 | 5 | A |
| LCR062R4 ^A | 66.0 x 33.0 x 97.0 | 104.0 | 6v./2.4Ah | 0.56 | Now 04204 (WP2.8-6P) | 8 | A |
| LCR063R2 ^A | 66.0 x 33.0 x 118.0 | 125.0 | 6v./3.8Ah | 0.73 | Now 04202 (WP3.8-6P) | 8 | A |
| LCR063R4 | 134.0 x 34.0 x 60.0 | 66.0 | 6v./3.4Ah | 0.62 | 04212 | 5 | A |
| LCR064R5P | 70.0 x 48.0 x 102.0 | 108.0 | 6v./4.5Ah | 0.72 | 04282 | 9 | A |
| LCR067R2P | 151.0 x 34.0 x 94.0 | 100.0 | 6v./7.2Ah | 1.26 | 04273 | 1 | A |
| LCR0612P | 151.0 x 50.0 x 94.0 | 100.0 | 6v./12.0Ah | 1.95 | 04232 | 1 | A |
| LCR121R3PG | 97.0 x 47.5 x 50.0 | 55.0 | 12v./1.3Ah | 0.59 | 04228 | 3 | A |
| LCR122R2PG | 177.0 x 34.0 x 60.0 | 66.0 | 12v./2.2Ah | 0.80 | 04200 | 1 | A |
| LCR123R4PG | 134.0 x 67.0 x 60.0 | 65.0 | 12v./3.4Ah | 1.20 | 04233 | 3 | A |
| LCR127R2PG | 151.0 x 64.5 x 94.0 | 100.0 | 12v./7.2Ah | 2.47 | 04274 | 4 | A |
| LCR127R2P1 | 151.0 x 64.5 x 94.0 | 100.0 | 12v./7.2Ah | 2.47 | 04224 | 4 | D |
| LCRA1212P1 | 151.0 x 98.0 x 94.0 | 102.0 | 12v./12.0Ah | 3.80 | 04211 | 4 | D |
| LCXD1217P | 181.0 x 76.0 x 167.0 | 167.0 | 12v./17.0Ah | 6.50 | 04235 | 2 | POST ^C |
| LCX1228AP | 165.0 x 125.0 x 175.0 | 179.5 | 12v./28.0Ah | 11.0 | 04220 | 2 | POST ^D |
| LCS122AG ^A | 201.0 x 25.0 x 61.0 | N/A | 12v./2.0Ah | 0.72 | Now 04201 (WP1220) | 10 | - |
| LCSA122R3AG ^{AB} | 182.0 x 24.0 x 62.0 | N/A | 12v./2.3Ah | 0.75 | Now 04231 (12v./2.3Ah) | 11 | - |
| UPRW1245P | 151.0 x 64.5 x 94.0 | 101.5 | 12v./45.0w | 2.60 | 04287 | 4 | A |
| UPRWA1232P1 | 151.0 x 51.0 x 94.0 | 100.0 | 12v./32.0w | 2.00 | 04292 | 4 | D |

^A Obsolete - an alternative battery (non Panasonic) is available. See page 26. Use part numbers to cross reference details.

^B VJ0180 Mating connector available for LCS type batteries (part no: 85174).

^C Terminal posts with nut & bolt connectors.

^D M5 threaded post.

PLEASE ENSURE YOU CAREFULLY CHECK BATTERY SPECIFICATIONS AND THE LAYOUT AND TERMINAL INFORMATION ON PAGE 29. ALL BATTERIES HIGHLIGHTED ABOVE ARE OF THE ABSORBENT GLASS MATT (AGM) TYPE.

Other Sealed Lead Acid Batteries

| MODEL | L x W x H (mm) | HEIGHT TO TERMINALS (mm) | VOLTAGE/CAPACITY | WEIGHT APPROX. (kg) | PART NO. | LAYOUT | TERMINAL |
|----------------------------|-----------------------|--------------------------|------------------|---------------------|---------------------|--------|--------------------------------|
| DRYFIT (SONNENSCHN) | | | | | | | |
| A204/1.0K | 34.9 x 42.5 x 50.5 | 51.7 | 4v./1.0Ah | 0.18 | 04207 | 13 | CONTACT |
| A506/1.2S | 97.3 x 25.5 x 51.0 | 54.9 | 6v./1.2Ah | 0.33 | 04208 | 14 | A |
| A506/3.5S | 134.3 x 34.8 x 60.5 | 64.4 | 6v./3.5Ah | 0.70 | 04215 | 15 | A |
| A208/2.5S | 133.5 x 36.5 x 63.0 | 66.6 | 8v./2.5Ah | 0.78 | 04255 | 16 | A |
| A508/3.5S | 178.5 x 34.1 x 60.5 | 64.4 | 8v./3.5Ah | 0.98 | 04242 | 17 | A |
| A512/2.0S (A212/1.8S) | 178.5 x 34.1 x 60.5 | 64.4 | 12v./2.0Ah | 1.00 | 04236 | 18 | A |
| A512/3.5S (A212/3.0S) | 134.5 x 66.8 x 60.5 | 64.4 | 12v./3.5Ah | 1.50 | 04223 | 19 | A |
| A512/6.5S | 151.7 x 65.5 x 94.5 | 98.4 | 12v./6.5Ah | 2.60 | 04267 | 20 | A |
| A512/10S | 152.0 x 98.0 x 94.5 | 98.4 | 12v./10.0Ah | 2.10 | 04248 | 21 | A |
| GF12033Y1 (A512C/38A) | 193.0 x 175.0 x 175.0 | 175.0 | 12v./38.0Ah | 15.00 | 04214 | 22 | POST ⁶ |
| EURO ENERGY | | | | | | | |
| 12v./2.3Ah | 182.0 x 24.0 x 62.0 | N/A | 12v./2.3Ah | 0.75 | 04231 | 11 | - |
| FIAMM | | | | | | | |
| FG10301 | 134.0 x 34.0 x 60.0 | 67.0 | 6v./3.0Ah | 0.70 | 04113 | 1 | A |
| FG20121 | 97.0 x 48.5 x 50.5 | 55.0 | 12v./1.2Ah | 0.62 | 04257 | 4 | A |
| FG20121A | 98.0 x 42.0 x 52.0 | 58.0 | 12v./1.2Ah | 0.55 | 04272 | 3 | A |
| FG20291 | 132.0 x 33.0 x 98.0 | 104.0 | 12v./2.9Ah | 1.17 | 04271 | 1 | A |
| GLOBAL YUASA | | | | | | | |
| ES4-12 ^A | 195.0 x 47.0 x 70.0 | 74.0 | 12v./4.0Ah | 1.70 | 04219 | 1 | A |
| HITACHI | | | | | | | |
| HV2.6-12 ^A | 178.0 x 34.0 x 60.0 | 66.0 | 12v./2.6Ah | 1.08 | Now 04293 (PX12026) | 4 | A |
| PBQ | | | | | | | |
| PBQ4.5-6 | 70.0 x 47.0 x 101.0 | 105.0 | 6v./4.5Ah | 0.85 | 04294 | 5 | A |
| PBQ1.3-12 | 98.0 x 48.0 x 58.0 | 58.0 | 12v./1.3Ah | 0.56 | 04254 | 4 | A |
| PBQ2.9-12 | 79.0 x 55.5 x 99.0 | 102.0 | 12v./2.9Ah | 1.17 | 04249 | 2 | A |
| PBQ4-12 | 195.0 x 47.0 x 70.0 | 74.0 | 12v./4.0Ah | 1.70 | 04219 | 1 | A |
| PORTALAC | | | | | | | |
| PX12026 | 178.0 x 34.0 x 60.0 | 66.0 | 12v./2.6Ah | 1.08 | 04293 | 4 | A |
| POWERSONIC | | | | | | | |
| PS-682 | 98.0 x 56.0 x 118.0 | 120.0 | 6v./8.0Ah | 1.52 | 04227 | 12 | A |
| PS-6120FP | 108.0 x 71.0 x 141.0 | 141.0 | 6v./12.0Ah | 2.40 | 04290 | 1 | D ^E /A ^F |
| PS-832 | 132.0 x 35.0 x 66.0 | 69.0 | 8v./3.2Ah | 0.86 | 04278 | 3 | A |
| PS-1208 | 96.0 x 25.0 x 62.0 | N/A | 12v./0.8Ah | 0.35 | 04276 | 6 | I |
| PS-1229 | 178.0 x 34.0 x 60.0 | 66.0 | 12v./2.9Ah | 1.08 | 04216 | 1 | A |
| PS-12330 | 197.0 x 132.0 x 165.0 | 180.0 | 12v./33.0Ah | 10.50 | 04259 | 1 | POST ⁶ |
| UNIPOWER | | | | | | | |
| WP2.8-6P | 66.0 x 33.0 x 97.0 | 104.0 | 6v./2.8Ah | 0.56 | 04204 | 8 | A |
| WP3.8-6P | 66.0 x 33.0 x 118.0 | 125.0 | 6v./3.8Ah | 0.73 | 04202 | 8 | A |
| WP1250 | 143.0 x 23.0 x 64.0 | N/A | 12v./2.0Ah | 0.65 | 04265 | 11 | - |
| WP1220 | 201.0 x 25.0 x 61.0 | N/A | 12v./2.0Ah | 0.72 | 04201 | 10 | - |
| WP2.6-12 ^A | 178.0 x 34.0 x 60.0 | 66.0 | 12v./2.6Ah | 1.08 | Now 04216 (PS-1229) | 1 | A |

^A Obsolete - an alternative battery is available. Use part numbers to cross reference details on this page.


^C Terminal posts with nut & bolt connector.

^E Positive terminal.

^F Negative terminal.

⁶ SAE Post.

DRYFIT/SONNENSCHN BATTERIES HIGHLIGHTED ABOVE ARE OF THE GEL VARIETY. OTHER SEALED LEAD ACID BATTERIES LISTED ABOVE ARE OF THE ABSORBENT GLASS MATT (AGM) TYPE.




Emergency Lighting Batteries

TYPE 1




TYPE 3



TYPE 7



TYPE 2




TYPE 4



TYPE 8



TYPE 5



| TYPE | DESCRIPTION | VOLTAGE/CAPACITY | PART NO. |
|------|--|------------------|----------|
| 1 | 3 Cell Stick With Solder Tags | 3.6v./4.0Ah | 01723 |
| 1 | 2 Cell Stick With Solder Tags | 2.4v./4.0Ah | 01722 |
| 2 | 3 Cell Stick With End Caps, Leads & Connector ^{A D} | 3.6v./4.0Ah | EE170030 |
| 2 | 2 x 3 Cell Stick With Single End Caps, Link Wire & Leads ^C | 7.2v./4.0Ah | EE170132 |
| 3 | 2 x 2 Cell Stick With Double Ends Caps, Leads & Connector ^{A D} | 4.8v./4.0Ah | EE170253 |
| 1 | 3 Cell Stick With Leads ^B | 3.6v./4.0Ah | EE170059 |
| 1 | 4 Cell Stick With Solder Tags | 4.8v./4.5Ah | 01724 |
| 2 | 4 Cell Stock With Single End Caps, Leads & Connector ^{A D} | 4.8v./4.0Ah | EE170038 |
| 3 | 2 x 3 Cell Stick With Double End Caps, Leads & Connector ^{A D} | 7.2v./4.0Ah | EE170262 |
| 1 | 2 Cell Stick With Fly Leads ^B | 2.4v./4.0Ah | EE170101 |
| 2 | 2 Cell Stick With Single End Caps, Leads & Connector ^{A D} | 2.4v./4.0Ah | EE170032 |
| 5 | 3 Cell Row with Backplate, Leads & Connector ^{A D} | 3.6v./4.0Ah | EE170041 |
| 4 | 2 x 2 Cell Stick With Double Ends Caps, Leads & Connector ^{A D} | 4.8v./4.0Ah | EE170285 |
| 5 | 2 Cell Row With Backplate, Leads & Connector ^{A D} | 2.4v./4.0Ah | EE170204 |
| 1 | 5 Cell Stick With Solder Tags | 6v./4.5Ah | 01725 |
| 2 | 5 Cell Stick With Single End Caps, Leads & Connector ^{A D} | 6v./4.0Ah | EE170008 |
| 1 | 4 Cell Stick With Leads ^B | 4.8v./4.0Ah | EE170105 |
| 5 | 4 Cell Row with Backplate, Leads & Connector ^{A D} | 4.8v./4.0Ah | EE170048 |
| 7 | 3 Cell Stick With Solder Tags ^E | 3.6v./1.2Ah | EE120037 |
| 2 | 6 Cell Stick With Single End Caps, Leads & Connector ^{A D} | 7.2v./4.0Ah | EE170053 |
| 8 | 10 Cell Block With Top, Bottom Cap & Solder Tags ^E | 12v./1.2Ah | EE120313 |

A 200mm Leads

B 300mm Leads

C 350mm Leads

D Standard Emergency Lighting Connector

E C Cells

All Batteries Used are D Cells Unless Indicated

OTHER EMERGENCY LIGHTING BATTERY COMBINATIONS ARE AVAILABLE. PLEASE CALL OUR CUSTOMER SERVICE TEAM ON (0116) 234 0567 FOR FURTHER INFORMATION.

Battery Chargers



UK PLUG TOP CHARGER
PART Nos: 04984, 04968, 04985



TWIN BENCH TOP BATTERY CHARGER
PART No: 04952

Chargers for Sealed Lead Acid Batteries

Multi-stage charging sequence approved to give maximum service life in all SLA batteries

- Orange LED – bulk charge
- Green LED – ready for use

| MODEL | STYLE | VOLTAGE | OUTPUT | BATTERY CAPACITIES | PART NO. |
|--------------|-------------|---------|--------|--------------------|----------|
| 9641-6/2.5 | UK Plug Top | 6v. | 2.5A | 5.0 - 24.0Ah | 04984 |
| 9641-12/1.0 | UK Plug Top | 12v. | 1.0A | 2.3 - 5.0Ah | 04968 |
| 9641-12/2.7 | UK Plug Top | 12v. | 2.7A | 5.0 - 24.0Ah | 04985 |
| 2040-12/4.0 | Bench Top | 12v. | 4.0A | 12.0 - 48.0Ah | 04993 |
| 2040-24/2.0 | Bench Top | 24v. | 2.0A | 5.0 - 12.0Ah | 04987 |
| 9740-12/10.0 | Bench Top | 12v. | 10.0A | 38.0 - 100.0Ah | 04990 |
| 9740-24/5.0 | Bench Top | 24v. | 5.0A | 12.0 - 48.0Ah | 04970 |

Chargers for Panasonic LCS/T Batteries (SLA Batteries)

| MODEL | STYLE | VOLTAGE | BATTERY CAPACITIES | PART NO. |
|--|-----------------|---------|--------------------|----------|
| LCS122AG/LCSD122PG/LCSA122R3AG/LCTA121R9PG | Single Plug Top | 12v. | 2.0 - 2.3Ah | EE460034 |
| LCSA122R3AG | Twin Bench Top | 12v. | 2.0 - 2.3Ah | 04952 |


Checkmate III Battery Capacity Tester

For Ni-Cd, Ni-MH, VRLA, Cyclon and Li-ion batteries

| MODEL | PART NO. |
|---|----------|
| Checkmate III Test Unit | 04957 |
| Checkmate III Power Supply | 04956 |
| Checkmate III Printer | 04958 |
| Checkmate III Printer Paper (10 rolls of thermal paper) | 04960 |
| Recalibration of Checkmate III test unit, checking safety of power supply and calibration certificate | 70004 |
| Recalibration of Checkmate III test unit only and calibration certificate | 70005 |
| Upgrade Checkmate II to Checkmate III | 70006 |

- Checkmate III allows data to be exported to PC
- Checkmate II upgrade available

PLEASE NOTE THAT OTHER CHARGERS ARE AVAILABLE, CALL OUR CUSTOMER SERVICE TEAM ON (0116) 234 0567 FOR FURTHER INFORMATION.



Battery Terminal & Layout Information

Terminals

| <p>Faston tab: 187 A</p> <table border="1"> <thead> <tr> <th colspan="2">INCH = MM</th> </tr> </thead> <tbody> <tr><td>.250</td><td>6.35</td></tr> <tr><td>.185</td><td>4.70</td></tr> <tr><td>.124</td><td>3.15</td></tr> <tr><td>.098</td><td>2.50</td></tr> <tr><td>.059</td><td>1.50</td></tr> <tr><td>.031</td><td>0.80</td></tr> <tr><td>.020</td><td>0.50</td></tr> <tr><td>.004</td><td>0.10</td></tr> </tbody> </table> | INCH = MM | | .250 | 6.35 | .185 | 4.70 | .124 | 3.15 | .098 | 2.50 | .059 | 1.50 | .031 | 0.80 | .020 | 0.50 | .004 | 0.10 | <p>Faston tab: 187 B</p> <table border="1"> <thead> <tr> <th colspan="2">INCH = MM</th> </tr> </thead> <tbody> <tr><td>0.472</td><td>12.00</td></tr> <tr><td>0.250</td><td>6.35</td></tr> <tr><td>0.236</td><td>6.00</td></tr> <tr><td>0.185</td><td>4.70</td></tr> <tr><td>0.079</td><td>2.00</td></tr> <tr><td>0.020</td><td>0.50</td></tr> </tbody> </table> | INCH = MM | | 0.472 | 12.00 | 0.250 | 6.35 | 0.236 | 6.00 | 0.185 | 4.70 | 0.079 | 2.00 | 0.020 | 0.50 | <p>Faston tab: 250 D</p> <table border="1"> <thead> <tr> <th colspan="2">INCH = MM</th> </tr> </thead> <tbody> <tr><td>.310</td><td>7.90</td></tr> <tr><td>.250</td><td>6.35</td></tr> <tr><td>.16</td><td>4.0</td></tr> <tr><td>.031</td><td>0.8</td></tr> <tr><td>.020</td><td>0.5</td></tr> </tbody> </table> | INCH = MM | | .310 | 7.90 | .250 | 6.35 | .16 | 4.0 | .031 | 0.8 | .020 | 0.5 |
|---|-----------|----|------|-------|------|-------|------|-------|------|------|-------|------|--|-----------|------|------|------|------|--|-----------|-----|-------|-------|--|----------|-------|------|-------|------|-------|------|-------|------|---|-----------|----|------|------|------|------|-----|-----|------|-----|------|-----|
| INCH = MM | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .250 | 6.35 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .185 | 4.70 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .124 | 3.15 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .098 | 2.50 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .059 | 1.50 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .031 | 0.80 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .020 | 0.50 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .004 | 0.10 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| INCH = MM | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.472 | 12.00 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.250 | 6.35 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.236 | 6.00 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.185 | 4.70 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.079 | 2.00 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.020 | 0.50 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| INCH = MM | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .310 | 7.90 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .250 | 6.35 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .16 | 4.0 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .031 | 0.8 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .020 | 0.5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <p>Bolt fastened terminal E</p> <table border="1"> <thead> <tr> <th colspan="2">INCH = MM</th> </tr> </thead> <tbody> <tr><td>.472</td><td>12.00</td></tr> <tr><td>.453</td><td>11.50</td></tr> <tr><td>.433</td><td>11.00</td></tr> <tr><td>.216</td><td>5.50</td></tr> <tr><td>.079</td><td>2.00</td></tr> </tbody> </table> | INCH = MM | | .472 | 12.00 | .453 | 11.50 | .433 | 11.00 | .216 | 5.50 | .079 | 2.00 | <p>Bolt fastened terminal F</p> <table border="1"> <thead> <tr> <th colspan="2">INCH = MM</th> </tr> </thead> <tbody> <tr><td>.59</td><td>15.0</td></tr> <tr><td>.216</td><td>5.5</td></tr> <tr><td>.295</td><td>7.5</td></tr> <tr><td>.197</td><td>5.0</td></tr> </tbody> </table> | INCH = MM | | .59 | 15.0 | .216 | 5.5 | .295 | 7.5 | .197 | 5.0 | <p>Internal thread terminal</p> <table border="1"> <thead> <tr> <th>Terminal</th> <th>x</th> <th>y</th> </tr> </thead> <tbody> <tr><td>J:</td><td>17mm</td><td>M5</td></tr> <tr><td>K:</td><td>20mm</td><td>M6</td></tr> <tr><td>L:</td><td>20mm</td><td>M8</td></tr> </tbody> </table> | Terminal | x | y | J: | 17mm | M5 | K: | 20mm | M6 | L: | 20mm | M8 | | | | | | | | | | |
| INCH = MM | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .472 | 12.00 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .453 | 11.50 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .433 | 11.00 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .216 | 5.50 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .079 | 2.00 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| INCH = MM | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .59 | 15.0 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .216 | 5.5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .295 | 7.5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| .197 | 5.0 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Terminal | x | y | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| J: | 17mm | M5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| K: | 20mm | M6 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| L: | 20mm | M8 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <p>Bolt fastened terminal G</p> <table border="1"> <thead> <tr> <th colspan="2">INCH = MM</th> </tr> </thead> <tbody> <tr><td>0.24</td><td>6</td></tr> <tr><td>0.27</td><td>7</td></tr> <tr><td>0.37</td><td>9.5</td></tr> <tr><td>0.71</td><td>18</td></tr> <tr><td>0.768</td><td>19.5</td></tr> </tbody> </table> | INCH = MM | | 0.24 | 6 | 0.27 | 7 | 0.37 | 9.5 | 0.71 | 18 | 0.768 | 19.5 | <p>JST No. VHR-2N I</p> <p>WIRE AWG #20 UL 1007</p> <p>JST No. VHR-2N (TERMINAL) JST SVH-21T-P1.1</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| INCH = MM | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.24 | 6 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.27 | 7 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.37 | 9.5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.71 | 18 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0.768 | 19.5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Layouts

1 2 3 4 5 6

7 8 9 10 11 12

Terminals/Layouts - Sonnenschein

13 14 15 16 17

18 19 20 21 22

Battery Care

Handling

- Avoid metal surfaces, metal tools, bracelets, jewellery, wire leads, etc. that may accidentally short-circuit the battery terminals.
- Do not drop or mechanically deform the battery as this may short-circuit the battery internally and cause fire.
- It is recommended that the battery be checked for mechanical damage before use.
- Damaged batteries may leak a corrosive fluid. If damaged, handle with protective goggles and gloves. Avoid contact with eyes and skin.
- Do not simultaneously handle both terminals of the battery. Under normal conditions, a battery voltage less than 40V is not considered dangerous. A voltage above 40V can be fatal.

Storage

- Store all batteries in a cool, dry environment (0-25°C preferred) away from flammable vapours and gases.
- Lithium Ion, Nickel Metal Hydride and Sealed Lead Acid batteries will self-discharge, but must not be allowed to become fully discharged. It is recommended that these batteries be recharged every six months.
- Nickel Cadmium batteries will slowly self-discharge with time, but may be stored either charged or discharged without damage.

Charging/Testing

- Battery charging and discharge testing must only be carried out by trained technicians following manufacturers instructions.
- Batteries will release hazardous materials and/or explode if overcharged, over-discharged or damaged.
- Pre-use checks should be carried out when batteries are used for critical devices.

Installation


- Each battery is designed for a specific use. If in doubt, check the specification number using either the current Euro Energy catalogue or by contacting our Customer Service team.
- Checks should be made to ensure that the correct battery is used in the relevant medical device.
- Batteries for use in critical devices should be checked and installed by a trained technician.
- Batteries should only be changed by personnel trained to operate the medical device being powered.

Shelf Life

- Rechargeable batteries supplied by Euro Energy are recommended to be installed within two years of the date of manufacture, subject to Sealed Lead Acid, Nickel Metal Hydride and Lithium Ion batteries being stored in a cool environment and fully recharged once every six months.
- It is recommended that a minimum of at least two batteries are kept in stock to ensure that a fully charged battery is ready for application at all times.
- Batteries have a finite life and therefore a log should be kept detailing the typical number of charge cycles.

Safe Disposal of Batteries

- Batteries must not be disposed with normal waste in shallow land fill as this contravenes the Waste Batteries & Accumulators Regulations 2009.(2009 No.890).
- Waste batteries for disposal must only be collected/transported by a carrier holding a Waste Management licence. Transportation by any other carrier contravenes the Special Waste Regulations 1996 and the Environmental Protection (Duty of Care Regulations) 1991. Offenders are liable to either or both a fine of a maximum of £20,000 and imprisonment of a maximum of five years.
- For battery waste enquires you may wish to contact G&P Batteries for advice on (0121) 568 3200.



Important Information

How to Order

This catalogue will help you to rapidly and easily select the product you need. There are a variety of ways in which to order. Please select one of the methods below when placing your order:

Telephone: (0116) 234 0567

Fax: (0116) 235 4183

Email: orders@euroenergy.co.uk

Post: Euro Energy Resources Limited, 10 Barrington Business Park, Leycroft Road, Leicester, LE4 1ET

Note that we can accept the following credit cards – Visa, Mastercard, Delta, Solo, and Switch/Maestro.

Always help us to provide the best service possible by quoting your account number, where relevant, when placing your order.

Please contact us for carriage charges.

VAT Exempt Supplies

All requests for the VAT exempt supply of replacement medical batteries must be supported by an appropriate exemption certificate. When ordering, please ensure that VAT exemption is specified and that the supporting certificate is forwarded to us at the earliest opportunity. Retrospective claims for VAT exemption may result in an additional administration charge being imposed. No VAT exempt goods can be despatched unless and until a valid, signed certificate has been received.

Technical Information

Euro Energy battery packs are assembled by a team of skilled technicians to BS EN ISO 9001: 2008 quality standards. Most replacement batteries carry the **CE** mark of conformity. Our replacement medical batteries are designed to meet or exceed the specification of the original equipment manufacturer.

Our in-house technical team are able to offer full customer support and discuss specifications with you and provide further information. We also have access to technical resources of the leading battery manufacturers.

Refurbishment

Euro Energy can only accept refurbishments where components and cells are deemed to be of a certain standard, and has the right to decline orders where these standards are not met.

Performance Guarantee

Our products carry a one year warranty, provided they are used and maintained in conjunction with the correct equipment, and in accordance with given instructions. Any products returned will be subject to rigorous testing. This is undertaken in order to maintain our high standards. Test data will be used to determine suitability of refund or replacement. Smart batteries, lithium products and Checkmate battery capacity testers, power supplies and printers carry a one year warranty from the date of purchase. Original Equipment Manufacturer (OEM) medical batteries carry standard manufacturer warranties.

Important Information Continued

Euro Energy Battery Labels

Every battery purchased, excluding smaller primary types, carries a silver Euro Energy label containing important information and a unique reference number for complete traceability. It is essential that this label is not removed from the product as this may render the product warranty invalid.

Spec Number: Use this code to reorder or enquire about the battery

Capacity: This is given in Ah or mAh

Nominal Voltage: For your information

Manuf. Date: The month and year of manufacture/despatch

LOT: A unique reference for traceability

 This logo indicates a product conforming to European standards

 Battery should be handled with care

  Separate Waste - Landfill not permitted



Index

Replacement Batteries

| | | | | | |
|--|----|--------------------------------|----|--|----|
| 3M | 02 | FUKUDA DENSHI | 08 | PETERS | 15 |
| ABBOTT | 02 | GE HEALTHCARE | 08 | PHILIPS | 15 |
| ACTAMED | 02 | GOULD | 08 | PHYSIO CONTROL | 16 |
| AEROGEN | 02 | GRASEBY | 09 | PICKER INTERNATIONAL | 16 |
| AEROSOL PRODUCTS/MG ELECTRIC | 02 | GRAYSTAD | 09 | POLYSTAN NORTH AMERICA | 16 |
| AGILENT TECHNOLOGY | 02 | GREINER | 09 | PPG | 16 |
| AIR SHIELDS | 02 | HAEMOCUE | 09 | PROPAQ/PROTOCOL | 16 |
| ALARIS | 02 | HEMOCHRON | 09 | PULMONETICS | 17 |
| ALCON | 02 | HEINE | 09 | PURITAN BENNETT | 17 |
| ALOKA | 02 | HELLIGE GMBH | 09 | QUEST MEDICAL | 17 |
| AMEDA | 03 | HEWLETT PACKARD | 09 | RADIOMETER | 17 |
| AND | 03 | HME | 10 | REFLOLUX | 17 |
| ARCOMED | 03 | HONEYWELL | 10 | REMEMI | 17 |
| ARJO | 03 | HORWELL | 10 | RESPIRONICS | 17 |
| ASTROMED | 03 | HOSPAL | 10 | RIGEL | 17 |
| ATMOS | 03 | HUNTLEIGH NESBITT EVANS | 10 | S&W VICKERS | 17 |
| B.T. | 03 | IMED | 10 | SCALE-TRONIX INC. | 17 |
| BARD | 03 | IMEX MEDICAL SYSTEMS | 11 | SCHILLER | 17 |
| BAXTER | 03 | IMPACT MEDICAL | 11 | SECA | 18 |
| BCI INTERNATIONAL | 04 | INFRASONICS | 11 | SEWARD MEDICAL | 18 |
| BOEHRINGER | 04 | INTERNATIONAL TECHNADYNE CORP. | 11 | SHERWOOD MED./AMERICAN HOSP. SUPPLIES | 18 |
| BRAUN | 04 | INVACARE | 11 | SEIMED | 18 |
| BREAS MEDICAL | 04 | IN-VIVO RESEARCH INC. | 11 | SIEMENS/BURDICK/SPACELABS/IVAC | 18 |
| BURDICK/SIEMENS/SPACELABS | 04 | IVAC/ALARIS | 11 | SIMED | 18 |
| CADD | 04 | KEELER OPTICAL | 11 | SONICAID MEDICAL | 18 |
| CAMBRIDGE HEART INC. | 04 | KONTRON | 11 | SPACELABS MEDICAL | 18 |
| CAMBRIDGE INST. | 04 | LAERDAL | 12 | SUZUKEN KENZ | 19 |
| CAMINO LABORATORIES | 04 | LIBERATOR | 12 | T.O.P. | 19 |
| CARDIAC SCIENCE | 04 | LIFECARE | 12 | TERUMO | 19 |
| CARDINAL HEALTH | 04 | LIKO | 12 | TOBY CHURCHILL | 19 |
| CARDIET | 04 | M D E | 12 | TORONTO | 19 |
| CARDIO RAPID | 05 | MAQUET | 12 | TRAVENOL | 19 |
| CARDIOLINE | 05 | MARDIS | 12 | TRIMED CORP. | 19 |
| CARDISUNNY | 05 | MARQUETTE | 12 | VERATHON MEDICAL | 19 |
| CAREFUSION | 05 | MASIMO | 12 | VICKERS | 19 |
| CME | 05 | MEDELA A.G. | 12 | VIOMEDIX | 19 |
| COLIN MEDICAL INSTRUMENTS | 05 | MEDICAL DIAGNOSTIC EQUIPMENT | 13 | VITALOG MONITORING | 19 |
| COROMETRICS MEDICAL SYSTEMS | 05 | MEDIX INFUSION SYSTEMS | 13 | VITALOGRAPH | 20 |
| COW & GATE | 05 | MEDTRONIC | 13 | WEINMANN | 20 |
| CRITICARE SYSTEM INC. | 06 | MEG | 13 | WELCH ALLYN | 20 |
| CRITIKON/JOHNSON & JOHNSON/GE HEALTHCARE | 06 | MERIAN | 13 | WELMED | 20 |
| DAMECA | 06 | MICRO MEDICAL LTD | 13 | ZOLL | 20 |
| DANBY | 06 | MINOLTA | 13 | | |
| DATASCOPE/MINDRAY | 06 | MOLIFT | 13 | | |
| DATEX | 06 | NELLCOR PURITAN BENNETT | 13 | | |
| DATREND SYSTEM | 07 | NESBITT & EVANS | 14 | | |
| DEVILBISS HEALTHCARE | 07 | NIDEK | 14 | Additional Products | |
| DIAGNOSTIC ULTRASOUND CORP. | 07 | NIHON KOHDEN | 14 | ALKALINE AND RECHARGEABLE BATTERIES/ SINGLE CELL CHARGERS/LITHIUM BATTERIES | 21 |
| DRAEGER | 07 | NONIN | 14 | UPS REPLACEMENT BATTERIES | 22 |
| DYLADE | 07 | NOVAMETRIX MEDICAL | 14 | CYCLON SEALED LEAD ACID BATTERIES | 23 |
| DYNATECH/NEVADA | 07 | OHMEDA MEDICAL PRODUCTS | 14 | YUASA SEALED LEAD ACID BATTERIES | 24 |
| EAST RADCLIFFE | 07 | OSCAR BOSCAROL | 15 | PANASONIC SEALED LEAD ACID BATTERIES | 25 |
| ELMED | 07 | OXFORD | 15 | OTHER SEALED LEAD ACID BATTERIES | 26 |
| ENGSTROM | 07 | PACETECH | 15 | EMERGENCY LIGHTING BATTERIES | 27 |
| ESCHMANN | 07 | PALCO LABORATORIES | 15 | SEALED LEAD ACID BATTERY CHARGERS/ CHECKMATE III BATTERY CAPACITY TESTER | 28 |
| FERNO | 08 | PARI | 15 | BATTERY TERMINAL & LAYOUT INFORMATION | 29 |
| FIBRE OPTIC MEDICAL | 08 | PARKS ELECTRONICS LABS | 15 | BATTERY CARE | 30 |
| FLAEM NUOVA | 08 | PEGASUS | 15 | HOW TO ORDER AND IMPORTANT INFORMATION | 31 |
| FRESENIUS | 08 | PENLON BROMSGROVE | 15 | | |



Euro Energy Resources Limited
10 Barrington Business Park
Leycroft Road, Leicester LE4 1ET

Tel: (0116) 234 0567
Fax: (0116) 235 4183

info@euroenergy.co.uk
www.euroenergy.co.uk

