

Labotron

The handy bench-scale orbital shaker



- ▶ Flexible
- ▶ Reliable
- ▶ Robust
- ▶ Corrosion resistant



www.infors-ht.com

May be subject to technical amendments

INFORS HT

Reliable shaking of small loads

► Flexible

Labotron is at home on the laboratory bench or in an incubator.

► Reliable

The sealed-for-life bearing lubricated with a high-temperature grease ensures a long and productive working life whatever the environment.

► Robust

The unit has a strong housing made from composite material. The shaking table is made from a solid plastic. This makes the whole easy to clean and resistant to abrasion and a wide range of chemicals.

► Corrosion resistant

The composite housing is corrosion resistant and ensures easy cleaning.



Technical Data

Shaker Speed	40–400 RPM
Throw	∅ 12,5 mm
Permissible	
Ambient T°	4°C up to 60°C
Ambient humidity	up to max. of 80° rF
Weight without table	Approx. 9 kg
Dimensions W x D x H	39 x 32 x 9 cm
Power supply	115/230V ± 10%, 50/60Hz

Loading Capacity

For operation using accessories without the need for the special screws, a table can be fitted.	Size of flask	Quantity
This is of type ATW 300 X 370 mm and can be universal, fixed or with adhesive matting.	50	35
	100	23
	250	12
	500	8

Infors AG
Headoffice, Switzerland

Rittergasse 27
CH-4103 Bottmingen
T +41 (0)61 425 77 00
F +41 (0)61 425 77 01
headoffice@infors-ht.com

Infors GmbH

Dachauer Str. 6
D-85254 Einsbach
T +49 (0)8135/83 33
F +49 (0)8135/83 20
infors.de@infors-ht.com

Infors UK Ltd

The Courtyard Business
Centre, Dovers Farm
Lonesome Lane, Reigate
Surrey, UK-RH2 7QT
T +44 (0)1737 22 31 00
F +44 (0)1737 24 72 13
infors.uk@infors-ht.com

Infors Sarl

6, rue Marcel Paul
FR-91742 Massy Cedex
T +33 (0)169 30 95 04
F +33 (0)169 30 95 05
infors.fr@infors-ht.com

Infors Benelux BV

Logistiekweg 20-02
7007 CJ Doetinchem
P.O. Box 613
NL-7000 AP Doetinchem
T + 31(0)314-36 44 50
F + 31(0)314-37 80 76
infors.bnl@infors-ht.com

Infors Canada

8350 rue Bombardier
Anjou, Quebec
Canada H1J 1A6
T + 1 514 352 50 95
F + 1 514 352 56 10
infors.ca@infors-ht.com

Infors Bio-Technology
(Beijing) Co., Ltd.

Room 505C, Building 106
Lize Zhongyuan
Wangjing New Industrial Zone
Chaoyang District, Beijing
100102 P.R. of China
T +86 10 51652068
F +86 10 64390585
info@infors-ht.com.cn