

CARLTON *Professional*



Beauty

Classic, Portable & Products

Carlton Professional Collection

We are pleased to present to you a well co-ordinated range of equipment, portable equipment and products designed by experts, built by Specialists and used by true Professionals worldwide.


Carlton Professional equipment is renowned for its high quality, durability and exceptional value. Used by Spas and Beauty Therapists for many years it offers service you have the right to expect.

We, at the Carlton Group, take pride in having been regularly chosen by Local Authorities, Colleges of Further Education and by virtually all Beauty Schools requiring updated or additional equipment or furniture. We also pride ourselves on the service we deliver to our clients remembering that they always have a choice, so we look forward to assisting you with your professional requirements.

Considerable care has been taken to ensure an accurate description in this brochure, but we reserve the right to change the specification of any item without prior notice. Please also note that sizes and colours may vary. For more information, please visit our website www.thecarltongroup.co.uk or contact our sales office staff who will be delighted to help you.



Angela Barbagelata
Chairman



Charlie-Claude Barbagelata
Managing Director

Portable Couch

With an all aluminium frame we have ensured that this couch is light to carry, easy to assemble, yet remarkably rigid for all treatments including full body massage. The adjustable back support has a sturdy single bar locking system and an integral breathe hole complete with upholstered plug. The frame is protected with a white electro-static powder coating and the top is sumptuously upholstered and covered with white reinforced vinyl. Each leg is secured by a lockable stay. Comes complete with a heavy duty and durable nylon carry bag with an "easy to carry" shoulder strap.

CC 950

SIZE: 1830 x 610 x 730mm
WEIGHT: 17.8kg



Magnifying Lamp

A professional folding magnifying lamp having a 3 diopter lens surrounded by a 22 watt circline bulb transformed to reduce heat output and produce a cool white light. For client comfort an anti-glare shield of translucent plastic is fitted around the lamp.

Supplied complete with a clamp and extension arm designed to fit directly on to the frame of your portable couch.

White, Chrome as standard.

Other colours special order.

Also available with a 5 diopter lens

REACH: 1200mm.



CC 100/LFAC

Portable Trolley

Two tier trolley. An attractive and practical trolley with large easy to steer castors.

Ideal for the mobile therapist or as an additional trolley for salon use. Unfolds, ready for use in seconds. White metal frame work and shelves.

Can accept up to 30kg.

CC 922

OPEN: 680 x 400 x 720mm
FOLDED: 680 x 90 x 980mm
WEIGHT: 6.2kg



CARLTON
Professional

Portable Epilation Portable Faradic



A solid state epilation unit utilising the latest technological advances in fully transistorised printed circuitry. After intensive research and development a system has been devised to operate at a specific frequency which will ensure excellent results with minimum trauma to the client.

The compact, attractive vanity case has spring clips fitted within the lid recess to store the switched needle holder and Swiss 'H' forceps that are provided with the unit. The front compartment has been constructed to hold several boxes of disposable needles together with the mains lead and moulded plug.

The fascia panel is fitted with a mains ON/OFF switch, a shock resistant fuse carrier, an outlet socket to connect the needle holder and an intensity control coupled to a visual indicator lamp.

CC 2345/E

SIZE: 350 x 250 x 130mm
WEIGHT: 3.2kg



This portable muscle exerciser uses the most modern transistorised circuitry with four outlets feeding eight conductive rubber electrodes through specially manufactured flat connecting leads. Each of these pairs of 3" electrodes has its own individual intensity control and operates totally independently.

To simplify operation the faradic muscle exerciser has a suitable frequency and pulse width pre-set in the factory during production.

In accordance with the current safety recommendations and to ensure that your client experiences no discomfort this unit incorporates a reset button which prevents its use unless all intensity controls have been set to zero prior to the commencement of a treatment. A facial electrode is supplied with the unit.

CC 2345/F

SIZE: 350 x 250 x 130mm
WEIGHT: 4.1kg

Portable High Frequency



This portable electronic high frequency unit is ideal for the mobile therapist, the vanity case having been constructed with spring clips in the recess of the lid to hold and protect the six electrodes that are supplied.

The sophisticated high frequency current changes 200,000 times per second and the electronic circuitry safety intensifies the voltage up to 150,000 volts with only minimal amperage. Because of this rapid oscillation there are no muscular contractions and the effect produced is either stimulating or soothing dependant on the chosen method of application.

CARLTON high frequency treatment is of tremendous benefit, its primary action being thermal, to produce heat. The current aids the deeper penetration of products and has a marked germicidal action on the skin.

The high frequency handle is of particular importance. The CARLTON handle is manufactured from material with a very high level of insulation to stop the current being transferred through to the therapist during a treatment.

CC 2345/HF

SIZE: 350 x 250 x 130mm
WEIGHT: 3.1kg

Portable HF Galvanic



The most popular combination available in a portable vanity case. The two functions of high frequency and galvanic treatment are contained in one compact unit that is supplied complete with full accessories for facial therapy treatment.

The MINIDERM has an electronic quartz minute timer to assist the therapist in planning her facial treatment programme. Each of the functions is controlled individually, and the fascia panel incorporates a liquid crystal display to monitor the output of galvanic current together with a polarity changeover switch.

Each unit is supplied complete with four glass high frequency electrodes, a galvanic desincrustor, mini-ball and indifferent electrode with sponge cover. To complete the accessory set is a pair of stainless steel galvanic rollers with nylon bushes and a separate link lead to connect them to the output from the unit.

To avoid confusion, the polarity changeover switch relates to the polarity of the black lead feeding the accessories, for example if the polarity switch is set to negative then any accessory connected to the black lead is negative.

CC 2345/MD

SIZE: 350 x 250 x 130mm
WEIGHT: 4.2kg

Portable Ultra Visage



Ultra Visage is now available in a portable vanity case. This microcurrent unit combines both galvanic with facial contour lifting for the optimum result. It is a low frequency mains operated electronic unit using pre-set timed variable wave forms.

Ultra Visage offers facial and body contour lifting by re-programming the muscle fibres, reducing and smoothing the fine lines and wrinkles on the neck. Sagging breasts and upper arms will improve immensely along with skin tone and texture.

CC 2345/UV

1 Set of products
1 Packet disposable sponges
1 Video CD

SIZE: 350 X 250 X 130mm
WEIGHT: 3.6kg



ChromaSolo



Used in combination with other treatments or separately, the Colour Therapy treatment stimulates the skin's cells to facilitate change. The inclusion of four colours (Red, Green, Blue and Yellow) using LED lights and set at specific wavelengths offers a wide range of benefits. For example, Green cools and reduces redness, Blue reduces swelling and offers anti-bacterial effects for the treatment of acne, Red rejuvenates and stimulates the lymphatic system, whilst Yellow tones and smoothes the skin, as well as promoting the production of collagen.

Four pulsed modes are also included which allow the wavelengths to penetrate the skin at different levels.

A highly versatile treatment, Colour Therapy can be incorporated into any treatment from dermabrasion to a massage.

CC 2385

SIZE: 100 x 165 x 55mm
WEIGHT: 525g



Audio Sonic Massage

This hand held, intra sound vibrator produces mixed frequency sound waves which penetrate up to 6cm into the body tissue and by increasing circulation help the body heal affected areas, quickly bringing relief from pain. Used in the treatment of arthritis, rheumatism, inflammation, circulation disorders, bronchitis, migraine and for the relief of nervous tension. Supplied with a 3 pin transformer plug together with two applicator heads.

CC 208

SIZE: 220 x 85 x 45mm

WEIGHT: 0.5kg



Vibromator

The Vibromator produces vibrations on the oscillation principle; that is it simultaneously gives out both horizontal and vertical vibrations. This action gives the Vibromator unusual versatility yet easy application whilst the two speed motor, for soft or firm massage, makes it equally suitable for the Beauty or Physio-therapist.

Supplied in a smart vanity box complete with a smooth rubber pad, a bumpy rubber pad and a plastic hyperaemia applicator.

PM 50

SIZE: 355 x 130 x 125mm

WEIGHT: 2.2kg



Lymph Drainage

Research has shown the need for an auto-matic facial massager to aid lymph drainage. The Vibrolift can be used over gauze on the face, or alternatively with your fingers through the vibrating loop to achieve a soothing hand massage.

CC 383

SIZE: 280 x 45mm

WEIGHT: 400g



CARLTON
Professional

High Frequency Facial Galvanic



An electronically generated high frequency current, controlled by an infinitely variable potentiometer. This allows the therapist to vary the output from a gentle soothing level, to a highly stimulating current, as required. Classic High Frequency has a well insulated H.F. handle and is supplied complete with both large and small facial bulbs and a glass saturator for indirect application.

HF 11

SIZE: 235 x 125 x 65mm
WEIGHT: 2.1kg



A facial electrotherapy unit for desincrustation and iontophoresis treatments. Electronic constant current ensures safe operation at all times by limiting the maximum output to the milliamp level that has been set. This output is clearly displayed on a digital read-out monitor. To avoid confusion each of the leads is colour coded and polarity is changed by a single switch that indicates which lead is the negative pole. Classic Facial Galvanic is supplied complete with stainless steel accessories comprising an indifferent, a desincrustor mushroom and a ball electrode together with connecting leads. Roller electrodes and link leads are available separately.

DC 06

SIZE: 235 x 125 x 65mm
WEIGHT: 2.1kg

Epilation



Literally a classic, this short wave diathermy unit uses solid state technology to produce very smooth current for user satisfaction and client comfort. The unit has been constructed for use with either footswitch or finger switch needle holder which are available to order.

Classic Epilation is adjusted by very accurately calibrated intensity control on the fascia panel and the current flow is monitored by an orange indicator lamp. This unit is very lightweight and strongly recommended for the visiting practitioner.

RF 08

SIZE: 235 x 125 x 65mm
WEIGHT: 2.0kg

Faradic



For tightening and improving muscle tone this neat unit incorporates two 25mm facial electrodes that may be detached from their moulded rubber carrier and used independently spaced when required and 4 x 3" pads for body work. Model FE10 has the pulse width and frequency pre-set in the factory and operator controlled relaxation, contraction length and intensity options. Client protection is provided by an integral re-set safety circuit. Supplied complete with non-tangle leads.

FE 10

SIZE: 235 x 125 x 65mm
WEIGHT: 2.1kg

Skin Care Plus

A range of preparations formulated for the care of skin in combination with GALVANIC and HIGH FREQUENCY treatments.

GALVANIC

- ELECTRO-LOTION
- DESINCRUSTANT GEL
- IONTO-GELS
- IONTO-AMPOULES
- MOISTURISING LOTION

HIGH FREQUENCY

- OZO CREAM

CARLTON Professional technology and gentle SKINCARE PLUS preparations enable high frequency and galvanic treatments to maximise the skin's potential.

SKINCARE PLUS preparations are particularly recommended for use with ULTRADERM, LIFT PROFESSIONAL, FACIAL GALVANIC, HIGH FREQUENCY or ULTRAPRO COMPLEX MACHINE.



Microdermabrasion

A comprehensive range of products to be used with PEEL PROFESSIONAL and BABYPEEL microdermabrasion machines.

- PRE-PEEL LOTION
- POST-PEEL LOTION
- ENZYMATIC PEELING PASTE
- OXYGEN MOISTURISER
- EGF SERUM
- EGF CREAM
- RICH HYALURONIC SERUM
- OPC REGENERATING CREAM
- VITAMIN C CREAM GEL
- CELL PROTECTING CREAM
- SKIN REFINER
- POST ABRASION SOOTHING GEL



Eyelash Tint



REKOCIL system was designed with your client in mind so as to ensure ultimate client comfort. This is achieved by using low levels of ammonia which means the product is especially mild. REKOCIL is very lightly perfumed thus avoiding the sensitisation usually associated with higher concentrations of fragrance.

REKOCIL comes in Black, Blue Black, and Brown. Supplied in 20gm tubes. Recommended for use with REKOCIL Special Developer and Eye Shields.

Needles

CARLTON Professional needles are manufactured with the latest state-of-the-art equipment which has been installed in a controlled 'clean-shop' and is operated by highly trained technicians under the supervision of skilled engineers with many years of experience in the manufacture of precision surgical instruments. The needles are diamond drawn from austenitic stainless steel, the steel used for hypodermic needles and surgical instruments, then beautifully finished with a specially developed final polishing process. When the needle is inserted into the shank it is crimped evenly on all three sides to give maximum assurance that the two pieces will not separate.

Wax



HOT WAX - Depilatory Wax - Original formula with pure beeswax. Also available with Azulene for use on a more sensitive skin. Supplied in 1kg bags.

WARM HONEY - Depilatory Wax - For safe, gentle removal of unwanted hair. Supplied in tin 400ml.

CARTRIDGE ROLLER - Depilatory Wax - three types of disposable cartridge waxes ready for use. 100ml cartridge.

CLASSIC NATURAL Honey Wax - all skin types.
LIPO GREEN Titanium dioxide Wax - for sensitive areas. Supplied with a small roller head.
GOLDEN PEARL Blue Ocean Wax - with Titanium Dioxide and Mica for the more delicate dry skin.



Naturalmoor



Naturalmoor therapeutic treatment system has been designed by European therapists, and is used in millions of applications each year.

Naturalmoor is an organo-mineralic-herbaceous matter that contains beneficial ingredients, which is found in central Europe in the high moors of the Czech Republic.

Naturalmoor organic matter is packed into porous sachets, it is indirectly heated using a state-of-the-art stainless steel Water Heater enabling the temperature to increase to an optimum it is then applied to the client's body. The heat is evenly distributed releasing the active ingredients to penetrate deeply, thus allowing the muscles to fully relax prior to the massage session. Pores dilate allowing the penetration of the organic mineral and herbaceous content, and allow the sweat to carry its bacteria and toxins back to the sachet.

Quite simply the Naturalmoor system is clean, simple and effective designed to achieve a treatment that is therapeutic, soothing and relaxing, based on a 20 minute heat treatment prior to massage, or just pure luxury during a facial treatment. Adding revenue but no extra time to the treatment.

- COMPLETELY NATURAL
- NO KNOWN SIDE EFFECTS
- SIMPLE TO APPLY
- ODOURLESS
- PAINLESS
- CLEAN
- RELAXING, STIMULATING
- ECOLOGICALLY FRIENDLY
- BIODEGRADABLE
- AFFORDABLE COURSE OF TREATMENT



Sterilsafe™ ... Totally New Technology

Sterilsafe™ ... The Complete System

Active against a wide range of organisms,
bacteria, algae, fungi, viruses...
(Including MRSA, SARS, TB, HEPATITUS B,
HIV-1, E-COLI, LEGIONELLA AND LISTERIA),
Medical Standard



COLD STERILSAFE SOLUTION

This formula containing biocides and polymer sterilises in ten minutes. The 1 Litre bottle of concentrate gives 40 weeks supply. Active solution is 100% effective for fourteen days.

- No Aldehydes
- Non Odorous
- Anti Rust Formula
- Concentrate used 19 : 1
- Convenient 1 Litre sterilising tray
- Non expensive peace of mind
- Non Hazardous
- Non Staining
- User Friendly



**Sterilsafe™
Spray**

As a ready to use spray it is convenient and always to hand. With a deodorizing action too it keeps areas cleaner for longer, giving total peace of mind.

500ml spray bottle.



**Sterilsafe™
Hand Mousse**

A simple, easy to use, pump action mousse to protect you and your clients. Un-perfumed - No alcohol. Totally cleanses the hands and forms a nano-spic film for long lasting protection. Gives rapid, long lasting protection.

600ml refillable bottle with pump.



Sterilsafe™

This indispensable hand cleanser and protector is also available in a practical 1 Litre container with handle in order to easily refill the 600ml hand mousse bottle.



**Sterilsafe™
Surface Wipes**

These very strong textured wipes are designed for thorough, effective cleaning of contaminated yet delicate surfaces. Non scratch wipe with Byotrol's unique anti-microbial formula.

Fragrance free. Alcohol free. Approx: 75 wipes.



The Carlton Group



TRAINING: At the Carlton Group our professionally qualified staff have an ongoing commitment to offer training, support and advice for all our equipment and products as a foundation for future development.

TRAINING CENTRE: Our 3000 sq.ft fully fitted training centre within the "Carlton Tower" enables the company to offer the finest post-graduate training and Continual Professional Development days, as well as courses on specific pieces of equipment. The centre offers prestigious facilities for up to 10 therapists to receive theoretical and practical training on specialised equipment in a professional yet relaxed atmosphere. Colleges are invited to arrange visits for lecturers and students to refresh and update their skills and information.



LUXURY SHOWROOM: Our luxurious showroom facilities have been designed as a spa/salon with the client's total comfort in mind. The showroom is presented as a series of individually stylised treatment rooms, each reflecting a different atmosphere through colour scheme and type of treatment, to offer the client a realistic environment in which to view, try out and analyse the equipment in a hands-on experience. Appointments are reserved to allow privacy and undivided attention.



EXPORT: With over 60 years of experience, Carlton Professional now exports to over 50 countries. The reputation and service offered by Carlton is renowned throughout the world. Used by many overseas Beauty Schools, Carlton equipment is requested by new therapists as their choice of equipment when opening their own salons. Supporting international exhibitions and visiting distributors to enhance their product knowledge and training contributes to increase our clientele all over the world. The "Carlton Family" overseas grows every day. Nowhere is too far and no job is too large or small. The personal attention given is our strength. *Everyone is important ...*



DELIVERY: As a manufacturer, the majority of our items are ex-stock, meaning our delivery service is prompt, efficient and reliable. Over the years our experience has enabled us to offer our clientele a high level of service and it is always our aim to constantly maintain this by *regularly updating our training and systems. We never forget you have a choice ...*

WARRANTY: Because of the confidence we, at the Carlton Group, have in our brand, we are very proud to offer a **five** year guarantee on all Carlton Professional electrotherapy machines, and two years on vapourisers and furniture. Should you ever experience a problem, there is always someone to help and advise you. *Remember we are only a call away ...*

CARLTON
Professional



CARLTON
Professional

For more information please visit
www.thecarltongroup.co.uk

For all UK orders and enquiries
Phone: 01903 761100 • Fax: 01903 751111



Carlton House, Commerce Way, Lancing,
West Sussex BN15 8TA, UK

Phone: +44 (0) 1903 761100 • Fax: +44 (0) 1903 751111

E-mail: info@thecarltongroup.co.uk • Web: www.thecarltongroup.co.uk